



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Basic Breast Examination Trainer, Light

Item No. ALT40202
Weight
Make: Adam, Rouilly Limited, UK

[Read More](#)

SKU:

Categories: Clinical Skills Trainers



Product Description

Limbs & Things Basic Breast Examination Trainer provides a highly realistic learning platform for acquiring the fundamental skills to perform Clinical Breast Examinations (CBE) and self breast examinations.

Skills

- Clinical Breast Examination (CBE)
- Self Breast Examination
- Identification of anatomical landmarks
- Location and diagnosis of pathologies
- Professional to patient communication with Simulated Patient set-up

Features

- Basic Breast Examination Trainer provides a highly realistic learning platform for acquiring the fundamental skills to perform Clinical Breast Examinations (CBE) and self breast examinations
- Interchangeable pathologies can be placed in various locations on the internal supports to represent different complications and pathologies
- Supplied with a variety of pathologies:
 - Cyst
 - Fibrocystic Disease
 - Fibroadenoma
 - 3 differently sized Carcinomas (2cm, 3cm and 5cm in size)
- This breast examination simulator was developed in collaboration with leading clinicians from Torbay Hospital, CSRC University of Liverpool and the Clinical Practice Centre, St James's University Hospital
- Pathologies are easily changeable and can be placed in various predetermined locations
- Soft tissue breasts offer a realistic look and feel
- Simulated Patient: wearable breast component can be used with a simulated patient, strap fastenings can be adjusted to the wearer
- Range of pathologies can be moved within the breast model for a variety of diagnostic scenarios
- Ergonomic design for Simulated Patient exercises, including prolonged use in OSCEs and assessments
- Adjustable straps to fit a wide range of simulated patients

Includes

- Wearable Examination Breasts
- ALT40204 Breast Examination Pathologies
- ALT40205 Breast Examination Inserts (Pair)
- ALT40206 Breast Pathologies Supports (Pair)
- ALT40207 Breast Back Plates (Pair)

Cleaning & Safety

- Product skin surface is washable using a mild soap and water
- Allow to fully dry before storing
- This product is latex free
- Use a small amount of talcum powder in the recess of the breast tissue to stop the insert sticking to the inside of the model
- When used on a simulated patient, ensure the breast trainer is fastened correctly, with not knots or twists in the straps

Manufacturer

Limbs & Things