



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Suprapubic Insertion Unit (Ultrasound Guided), Light

Item No. ALT60852
Weight
Make: Adam, Rouilly Limited, UK

[Read More](#)

SKU:

Categories:Ultrasound



Product Description

The Limbs & Things Suprapubic Insertion Unit can be used independently or combined with another catheterisation trainer for a real to life simulation. This module allows for the practice of the surgical procedure of suprapubic catheterisation and is made of echolucent silicone for ultrasound guided procedures.

Skills

- Suprapubic catheter management on female anatomy
- Blind and ultrasound guided suprapubic catheter insertion

- Ultrasound bladder scanning on a full bladder

Features

- Ultrasound bladder scanning – allows for ultrasound identification of urine present in the bladder, reducing risks associated with CaUTI
- Ultrasound guided Suprapubic Catheterisation – allows for advanced procedure training in Urology
- Allows for the practice of Blind Suprapubic Catheterisation
- More robust design means longer lasting models, keeping on-going costs lower
- Ability to view catheter path when modules are used outside of the Pelvic Shell
- Improved Water System
- Softer, more realistic anatomy
- Anatomy:
 - Pubic bone
 - Rectus sheath
 - Anatomically accurate bladder
- Improved internal anatomy to allow for smoother insertion of catheter into the bladder
- New ultrasound bladder scanning aides in identification of urine in the bladder, reducing risks associated with CaUTI
- Ultrasound guided catheterization allows for more advanced Urology training
- Modules are interchangeable, and fit with the Male and Female Urinary Catheterisation models
- Ability to view the path of the catheter during procedures
- Non-drip valve
- Simulated Patient: can be used with a standardised patient to help improve communication skills

Includes

- ALT60857 Suprapubic Catheterisation Module, Light
- ALT60858 Suprapubic Insertion Bungs (x2)
- ALT60863 Water System

Cleaning & Safety

- Drain all water from the product
- Remove all catheters from the trainer
- Clean product surface with soft, damp cloth and mild detergent
- Allow trainer to dry completely before storing
- Components of the trainer are latex free
- Always deflate the balloon before removing a catheter

Manufacturer

Limbs & Things