

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile :	+91 88495 63724	Email:	sem@semtrainers.com
Mobile 1:	+91 98791 03905	Website:	www.semtrainers.com



TruCric

Item No.	1021893
Weight	1.682kg
MPN:	TCRIC1
Read More	

SKU:

Categories: Airway Management Adult





Product Description

The TruCric is a realistic training model for cricothyrotomy, percutaneous tracheostomy and Air Jet ventilation procedures. Easily replaceable consumables make the TruCric Ideal for teaching

environments and facilitate repeated technique practice. Features include:

- Lightweight and durable training device designed for efficient procedure training
- Palpable cricoid landmarks, laryngeal cartilages and tracheal rings provide positive user feedback with easy identification of the sternal notch and the clavicle bone
- Correct technique verification with replaceable lung bag
- · Neck skin can be rotated on the model allowing up to 20 incisions
- Larynx insert allows 1 surgical cricothyroidotomy and 3 tracheostomy procedures
- All consumables are quick and easy to interchange, perfect for a busy teaching environment
- Optional subcutaneous fat tissue can be added to simulate difficult tracheostomy or cricothyroidotomy procedure

The TruCric will provide a solution to your training needs in:

- Needle and surgical cricothyroidotomy
- Air jet ventilation
- Percutaneous tracheostomy
- Seldinger technique

Product components:

- 1 Carrier case
- 1 USB pen drive user manual
- 1 Bottle of TruCorp® lubrication
- 1 Adult lung bag
- 10 Larynx inserts (1 attached and 9 spare)
- 2 Wrap around neck skins (1 attached and 1 spare)
- 1 Non-slip mat

	Deluxe Cricothyrotomy Simulator	TruCric	Cricotracheotomy Trainer	Cricothyrotomy Simulator
Item No.	<u>1017948</u>	<u>1021893</u>	<u>1005180</u>	<u>1020255</u>
Surgical cricothyrotomy	Yes	Yes	Yes	Yes
Needle cricothyrotomy	Yes	Yes	Yes	Yes
Insertion of endotracheal tube	Yes	Yes	Yes	Yes
Full use of devices such as the PCK and UPK	-	Yes	Yes	Yes
Seldinger technique	-	Yes	Yes	Yes
Anatomical Landmarks	Cricothyroid membrane	Cricothyroid membrane, laryngeal cartilages and tracheal rings	Cricothyroid membrane, chin and full-sized neck	Cricothyroid membrane, chin and full-sized neck
Review airway anatomy	-	-	Yes	Yes
Realistic soft skin	Yes	Yes	-	-

Confirm proper	Inflation of lung	Inflation of lung	no	Inflation of lung
placement of	indicates proper	indicates proper		indicates proper
endotracheal tube	placement	placement		placement
Simulation features	 Allows for an adult as well as pediatric cricothyrotomy Includes adult and 3- year-old tracheas 	• Neck skin can be rotated on the model allowing up to 20 incisions • Larynx insert allows 1 surgical cricothyroidotomy and 3 tracheostomy procedures	 Up to 36 incisions on each unit Up to 18 mini tracheostomies 	• Assessment of styled and obturator placement (replaceable trachea)

