



## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com



## Cricothyrotomy Simulator

**Item No.** LF00994  
**Weight** 6.7kg  
**Make:** Nasco Healthcare, USA

[Read More](#)

**SKU:**

**Categories:** Advanced Trauma Life Support (ATLS)

### Product Description

The newly designed **Life/form®** Cricothyrotomy Simulator has been developed for learning and practicing the techniques necessary to perform needle or surgical cricothyrotomy procedures. Paramedics, EMTs, combat medics, flight nurses, anesthesiologists, and other emergency medical personnel will have the opportunity to strengthen their ability and confidence to perform or assist in implementing surgical airways.

- Anatomically accurate landmarks aid in site training and allow for fast action.
- The hyperextended neck allows the user to determine the proper incision site. The trachea in the simulator is replaceable, as the airway passes completely through from top to bottom. This allows checking the stylet and obturator placement once the incision has been made.
- Complete with a chin and full-size neck, ties can be used to hold the obturator in a secure position. Inflation of the simulated lung verifies correct placement.

Includes the simulator with base, six replaceable neck skins, six adult trachea inserts (four rigid, two soft), six child trachea inserts (four rigid, two soft), two simulated lungs, instruction manual, and hard carry case.

	<b>DeLuxe Cricothyrotomy Simulator</b>	<b>TruCric</b>	<b>Cricotracheotomy Trainer</b>	<b>Cricothyrotomy Simulator</b>
<b>Item No.</b>	<a href="#">1017948</a>	<a href="#">1021893</a>	<a href="#">1005180</a>	<a href="#">1020255</a>
<b>Surgical cricothyrotomy</b>	Yes	Yes	Yes	Yes
<b>Needle cricothyrotomy</b>	Yes	Yes	Yes	Yes
<b>Insertion of endotracheal tube</b>	Yes	Yes	Yes	Yes
<b>Full use of devices such as the PCK and UPK</b>	-	Yes	Yes	Yes
<b>Seldinger technique</b>	-	Yes	Yes	Yes
<b>Anatomical Landmarks</b>	Cricothyroid membrane	Cricothyroid membrane, laryngeal cartilages and tracheal rings	Cricothyroid membrane, chin and full-sized neck	Cricothyroid membrane, chin and full-sized neck
<b>Review airway anatomy</b>	-	-	Yes	Yes
<b>Realistic soft skin</b>	Yes	Yes	-	-
<b>Confirm proper placement of endotracheal tube</b>	Inflation of lung indicates proper placement	Inflation of lung indicates proper placement	no	Inflation of lung indicates proper placement
<b>Simulation features</b>	<ul style="list-style-type: none"> <li>• Allows for an adult as well as pediatric cricothyrotomy</li> <li>• Includes adult and 3-year-old tracheas</li> </ul>	<ul style="list-style-type: none"> <li>• Neck skin can be rotated on the model allowing up to 20 incisions</li> <li>• Larynx insert allows 1 surgical cricothyroidotomy and 3 tracheostomy procedures</li> </ul>	<ul style="list-style-type: none"> <li>• Up to 36 incisions on each unit</li> <li>• Up to 18 mini tracheostomies</li> </ul>	<ul style="list-style-type: none"> <li>• Assessment of styled and obturator placement (replaceable trachea)</li> </ul>