



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Defibrillation Pad & Patient Adapter Physio Control Adapters for CRiSis Manikin

Item No. LF03656U

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:ALS Adult

Product Description

The reusable defibrillation pad adapters are designed to be used in place of expensive defib gel pads. They come complete with patient simulator adapters that will enable you to use the patient simulator you already have with the CRiSis™ Manikin, CRiSis™ Manikin Update Package, or Defibrillation Chest Skin. The training cables to use with the CRiSis™ manikin are available with unique cable ends for use with the following defibrillators: Zoll, Survivalink, Vivalink, and Medtronic Physio Control Quick combo. If your defibrillator uses snap-on patches we also have the adapters you need to eliminate the use of pads. Training cables are not recommended for use with non-rechargeable battery AED. This includes FirstSave, HeartStream, and Medtronic Physio control with non-rechargeable lithium batteries. This package will save you both time and money!