



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Simulator 12-Lead Arrhythmia with Manikin Overlay – Zoll, Large

Item No. 1020845
Brand Nasco
MPN: 101-12-000L

[Read More](#)

SKU:

Categories:BLS Adult

Product Description

Use on any large-sized manikin to change it into a 12-lead trainer. Pace and defibrillate directly on the

overlay system connected to the interactive 12-lead ECG simulator.

Training features:

- Four preset pacing capture levels: 70, 80, 90 or 100mA. Waveforms for pacing: Sinus Brady (two), 1st degree A-V block, 2nd degree type I A-V block, 2nd degree type II A-V block, 2nd degree type II A-V block with PVCs, and 3rd degree A-V block.
- Waveforms for defibrillator training: V. Fib, V. Tach (high rate), V. Tach (low rate), Torsade, A. Fib, A. Flutter, zPSVT, Sinus Tach, Sinus rhythm, Sinus rhythm with PVCs, and Asystole.
- ST segments and T wave of abnormalities: Anterior MI, Inferior MI, Antero-Septal ST elevation, Anterior ST depression, Lateral ST elevation, and Inferior ST elevation.