



## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com

## Advanced Child Airway Management Trainer

**Item No.** LF03762  
**Weight** 6kg  
**Make:** Nasco Healthcare, USA

[Read More](#)

**SKU:**

**Categories:** Airway Management Child



## Product Description

Present your students with the challenges they may face in the real world. Advanced Child Airway Management Trainers offer tongue swelling and laryngospasm. This advanced airway management trainer is perfect for practicing skills on pediatric patients and represents the size of an 5-year-old patient. Practice intubation, ventilation, suction, and jaw thrust techniques. Realistic anatomy and landmarks include teeth, tongue, oral and nasal pharynx, larynx, epiglottis, arytenoids, false and true vocal cords, trachea, lungs, esophagus, and stomach. The trainer allows you to practice oral, digital, and nasal intubation, as well as E.T. (Endotracheal Tube), E.O.A. (Esophageal Obturator Airway), P.T.L. (Pharyngeal Tracheal Lumen), L.M.A. (Laryngeal Mask Airway), E.G.T.A. (Esophageal Gastric Tube Airway), and Combitube® insertion. Separate lungs for auscultation. Inflatable stomach bladder indicates esophageal insufflation. With its slightly anterior position, swelling tongue, and vocal cords, the Advanced Child Airway Management Trainer is a great trainer for introductory as well as advanced

training. Includes training stand to mount head, pump spray lubricant, and hard carrying case.

	Head/torso	Age	Digital intubation	Tongue edema	Laryngospasm	Stomach rise	Suctioning	Manual carotid pulse	Defibrillation	Additional features
<b>NG Tube Insertion Trainer</b> <i>Item No: <a href="#">1017250</a></i>	head mounted on a board	1-year old	-	-	-	Yes	-	-	-	trainer for instruction of nasogastric tube insertion and care
<b>Infant Airway Simulator</b> <i>Item No: <a href="#">1018863</a></i>	full-body	1-year old	-	-	-	-	Yes	-	-	chest cavity containing realistic organs, nasal passage permits placement of NP tube
<b>Infant Airway Trainer</b> <i>Item No: <a href="#">1017954</a></i>	head mounted on a board	1-year old	-	-	-	Yes	Yes	-	-	accepts 4.0 mm endotracheal tubes
<b>Infant Intubation Head</b> <i>Item No: <a href="#">1017245</a></i>	head mounted on a board	1-year old	-	-	-	-	Yes	-	-	accepts 2.5 mm endotracheal tubes
<b>Advanced Infant Intubation Head</b> <i>Item No: <a href="#">1017236</a></i>	head mounted on a board	1-year old	-	Yes	-	Yes	Yes	-	-	Vocal cords are highlighted for easy viewing
<b>Airway Management Trainer</b> <i>Item No: <a href="#">1014201</a></i>	head mounted on a board	3-year old	-	-	-	-	Yes	-	-	accepts 4.0 mm endotracheal tubes

<b>Advanced Child Airway Management Trainer</b> <i>Item No:</i> <a href="#">1019787</a>	head	5-year old	-	Yes	Yes	-	Yes	Yes	-	can be attached to any of the Life/form® Child Airway Management Trainers or Child CRiSis Manikins
<b>Child Airway Management Trainer Head Item</b> <i>No:</i> <a href="#">1019788</a>	head	5-year old	-	-	-	Yes	Yes	Yes	-	can be attached to any of the Life/form® Child Airway Management Trainers or Child CRiSis Manikins
<b>Advanced Child Airway Management Trainer</b> <i>Item No:</i> <a href="#">1005698</a>	head mounted on a board	5-year old	Yes	Yes	Yes	Yes	Yes	Yes	-	also allows for E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and King System insertion
<b>Basic Child CPR/ Airway Management Torso"</b> <i>Item No:</i> <a href="#">1018865</a>	torso	5-year old	-	-	-	-	-	-	-	This torso is the foundation of the Child CRiSis™ system as arms, legs, and a defibrillation chest can be added later.

<b>Advanced Child CPR/Airway Management Torso</b> <i>Item No:</i> <a href="#">1018864</a>	torso	5-year old	Yes	Yes	Yes	Yes	Yes	Yes	Yes	also allows for E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and King System insertion
---	-------	------------	-----	-----	-----	-----	-----	-----	-----	---