



**SEM**

## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com



## Catheterization Simulator BASIC, male

Item No. 1020232 [P93B-M]  
Weight 2.3 kg  
Dimensions 36 x 11 x 26.5 cm  
Brand 3B Scientific

[Read More](#)

**SKU:**

**Categories:**Catheterization



## Product Description

With the male P93 BASIC catheterization trainer, male bladder catheterization procedures can be realistically demonstrated, practiced and assessed. The genital insert is placed into the holder and held in place with magnets. The material of the insert is soft and flexible: the foreskin is movable and the penis can be stretched. This way the student can practice all the necessary gestures for catheterization (e.g. disinfection). The material gives a realistic feeling when inserting and removing the catheter into

the bladder with realistic resistance and insertion depths.

**Direct learning:**

When the catheter has been correctly inserted, fluid runs out, just like with a real patient. The catheter can also be checked through the transparent bladder at any time.

**Various levels of difficulty:**

On the male insert, three different levels of narrowing of urethra can be set. In the narrowest setting, catheterization is not possible, leading to the conclusion that a suprapubic catheter is needed.

**All the benefits of the male P93 BASIC catheterization simulator at a glance:**

- Male catheterization with realistic resistance
- 3 levels of adjustable narrowing of urethra
- Soft and movable foreskin
- Liquid outflow if catheterization is successfully carried out
- The situation can be checked through the transparent bladder
- Stands securely on the table thanks to the non-slip feet
- Magnetic connectors for quick set-up and dismantling
- Easy to clean, and the inserts can be completely removed

**Delivery contents:**

- Stands with magnetic connections for genital inserts
- Catheterization set with transparent bladder
- Bottle and tube connectors

**Note:**

For catheterization, we recommend using catheters with only one or two drainage eyes opposite each other. Both silicone and PVC catheters can be used. We recommend the following catheters:

- Urethral narrowing “open” = 16 CH catheter or 14 CH balloon catheter
- Urethral narrowing “partial narrowing” = 14 CH catheter or 12 CH balloon catheter

**Catheters with a curved tip should not be used!****Tip:**

According to the guidelines for correct catheterization, even during the simulation, lubricating gel (not included) should be placed into the urethra to ensure that the catheter can be inserted smoothly.

Alternatively, some lubricating gel can be placed on the tip of the catheter.