

**SEM**

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Diabetic Injection Pad

Item No. WA18238U

Weight 0.179kg

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Intramuscular (I.m.) and Intradermal

Product Description

This pad that looks and feels like human skin is ideal for diabetic instruction and injection practice. May be worn on the body in specific locations (stomach, thigh, upper arm) to allow patients to inject fluid into the pad to practice giving themselves injections. Pad is thick enough to accept all insulin needles. For safety purposes, the 6-1/4" L x 4-1/4" W x 3/4" D oval pad has a hard plastic backing to prevent the needle from poking through the pad and into the wearer's skin. Has a 1" wide nylon strap with an adjustable plastic buckle.