



**SEM**

## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE


SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724

**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com

**Website:** www.semtrainers.com



## Intradermal Injection Simulator

Item No. 1005628 [W44097]

Weight 0.77 kg

Dimensions 33 x 28 x 13 cm

MPN: LF01008U

[Read More](#)

**SKU:**

**Categories:** Intramuscular (I.m.) and Intradermal

## Product Description

The simulator features a forearm from the wrist to just below the elbow. Vinyl skin provides realistic feel and appearance to ensure a realistic training experience. The simulator features eight sites for practicing intradermal injections. If fluid is properly injected, a characteristic skin welt will form. The welt is removed by withdrawing the fluid after practice. Each site is reusable by dozens of students.

Supplied with sealant, syringe and storage box.

	Intramuscular Injection Simulator - Upper Arm	Intramuscular Injection Simulator - Upper Leg	Buttock Injection Simulator	Intramuscular injection simulator	Intramuscular Injection Model of Buttocks	Intramuscular Injection Simulator	Intradermal Injection Simulator
<b>Item No.</b>	<a href="#">1009840</a>	<a href="#">1000511</a>	<a href="#">1000514</a>	<a href="#">1010008</a>	<a href="#">1005394</a>	<a href="#">1005586</a>	<a href="#">1005628</a>
<b>Body part</b>	upper arm	upper leg	buttock	muscle extract with skin	pelvis with buttocks	lower torso	forearm
<b>Palpable anatomical landmarks</b>	Yes	Yes	Yes	-	Yes	Yes	-
<b>Injection type</b>	intramuscular	intramuscular	intramuscular	intramuscular	intramuscular	intramuscular	intradermal
<b>Injection sites</b>	as indicated by a trainer	as indicated by a trainer	as indicated by a trainer	all the surface	as indicated by a trainer	as indicated by a trainer	eight
<b>Proper injection</b>	audiovisual feedback	audiovisual feedback	audiovisual feedback	-	green light will appear	-	skin welt will form
<b>Additional information</b>	training and test modes are available	training and test modes are available	training and test modes are available	will absorb fluid	special skin material does not show needle marks	simulated bony structure is embedded in the torso	liquid skin repair kit included (M-1013686)