



**SEM**

## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com

## Basic GERi™ Nursing Skills Elderly-Care Manikin

Item No. 1005620 [W44077]  
Weight 10.785 kg  
Dimensions 94 x 52 x 28 cm  
MPN: LF04040U

[Read More](#)

**SKU:**

**Categories:** Intramuscular (I.m.) and Intradermal



### Product Description

GERi™ has an elderly appearance with skin wrinkles and folds. Complete lifelike range of motion, realistic patient positioning, and non-pinching joints - moves like a person.

- Full-size adult manikin - measures 58".
- Lightweight - approximately 28 lbs.

- Overall female appearance
- Converts to male with removal of wig and attachment of male genitals

#### Visual Inspection:

- Normal and Cancerous Mole
- Stage 1 Sacral Ulcer
- Dilated and Constricted pupils
- Reddened Skin Folds

#### Patient Care Simulation:

- Bandaging and Wound Dressing
- Bed Baths
- Clothing Changes
- Denture Placement and Removal, Upper and Lower
- Ear Canal Irrigation, Otic Drops, and Hearing Aid Placement
- Eye Irrigation
- Finger and Toe Manipulation with Bandaging - Flexible and Individually Molded
- Hair Care Washing and Combing
- Intramuscular Injection Sites: Arm, Thigh, Buttock
- Oral and Nasal Hygiene: Lavage, Gavage, and Suctioning
- Ostomy Care: Ileostomy and Colostomy Tissue
- Maintenance and Appliance Application, Lavage, and Suctioning
- Patient Positioning
- Patient Transfer Techniques
- Pericare
- Tracheostomy Care, Lavage and Suctioning

Ideal for OBRA Required Training. Carry case not included.

|                                 | <b>Basic:</b><br>KERi™ <a href="#">1020187</a><br>GERi™ <a href="#">1005620</a> | <b>Complete:</b><br>KERi™ <a href="#">1013742</a><br>GERi™ <a href="#">1005597</a> | <b>Advanced:</b><br>KERi™ <a href="#">1017567</a><br>GERi™ <a href="#">1005606</a> | <b>Auscultation:</b><br>KERi™ <a href="#">1020147</a><br>GERi™ <a href="#">1020146</a> |
|---------------------------------|---------------------------------------------------------------------------------|------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|
| <b>Skin inspection and care</b> | Yes                                                                             | Yes                                                                                | Yes                                                                                | Yes                                                                                    |
| <b>Bandaging</b>                | Yes                                                                             | Yes                                                                                | Yes                                                                                | Yes                                                                                    |
| <b>Washing</b>                  | Yes                                                                             | Yes                                                                                | Yes                                                                                | Yes                                                                                    |
| <b>Denture Placement</b>        | Yes                                                                             | Yes                                                                                | Yes                                                                                | Yes                                                                                    |
| <b>Hearing aid placement</b>    | Yes                                                                             | Yes                                                                                | Yes                                                                                | Yes                                                                                    |

|                                            |     |     |     |     |
|--------------------------------------------|-----|-----|-----|-----|
| <b>Ear and eye irrigation</b>              | Yes | Yes | Yes | Yes |
| <b>Ostomy and tracheostomy care</b>        | Yes | Yes | Yes | Yes |
| <b>Positioning and transfer techniques</b> | Yes | Yes | Yes | Yes |
| <b>Pericare</b>                            | Yes | Yes | Yes | Yes |
| <b>Intramuscular injections</b>            | Yes | Yes | Yes | Yes |
| <b>Gastrostomy</b>                         | -   | Yes | Yes | Yes |
| <b>Pap smears and douching</b>             | -   | Yes | Yes | Yes |
| <b>Enema administration</b>                | -   | Yes | Yes | Yes |
| <b>Prostate exam</b>                       | -   | Yes | Yes | Yes |
| <b>Urinary catheterisation</b>             | -   | Yes | Yes | Yes |
| <b>Nasogastric tube placement</b>          | -   | Yes | Yes | Yes |
| <b>IV insertion</b>                        | -   | -   | Yes |     |
| <b>Blood pressure monitoring</b>           | -   | -   | Yes |     |
| <b>Recognition of heart sounds</b>         | -   | -   | -   | Yes |
| <b>Recognition of lung sounds</b>          | -   | -   | -   | Yes |