



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Buttock Injection Simulator

Item No. 1000514 [P57]

Weight 2.306kg

Dimensions 30 x 30 x 12 cm

Brand 3B Scientific

[Read More](#)

SKU:

Categories: Intramuscular (I.m.) and Intradermal



Product Description

This high-quality, strap-on intramuscular simulator is a lifelike model of a right buttock with all important anatomical landmarks for intramuscular (I.M.) injections: iliac crest, anterior superior iliac spine and greater trochanter. The integrated fine electronics provide extremely graphic feedback when training the intramuscular injection technique, correctly administered injections will produce audiovisual feedback. Additional training and control modes are available for checking results either immediately or

upon completion of a training session. This anatomically realistic injection simulator is supplied with a 21G/0.8 injection needle and a 5 ml syringe.

	Intramuscular Injection Simulator - Upper Arm	Intramuscular Injection Simulator - Upper Leg	Buttock Injection Simulator	Intramuscular injection simulator	Intramuscular Injection Model of Buttocks	Intramuscular Injection Simulator	Intradermal Injection Simulator
Item No.	1009840	1000511	1000514	1010008	1005394	1005586	1005628
Body part	upper arm	upper leg	buttock	muscle extract with skin	pelvis with buttocks	lower torso	forearm
Palpable anatomical landmarks	Yes	Yes	Yes	-	Yes	Yes	-
Injection type	intramuscular	intramuscular	intramuscular	intramuscular	intramuscular	intramuscular	intradermal
Injection sites	as indicated by a trainer	as indicated by a trainer	as indicated by a trainer	all the surface	as indicated by a trainer	as indicated by a trainer	eight
Proper injection	audiovisual feedback	audiovisual feedback	audiovisual feedback	-	green light will appear	-	skin welt will form
Additional information	training and test modes are available	training and test modes are available	training and test modes are available	will absorb fluid	special skin material does not show needle marks	simulated bony structure is embedded in the torso	liquid skin repair kit included (M-1013686)