



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



KERi™ Basic Manikin without Carry Case - Light

Item No. 1020187
Weight 11.2 kg
MPN: LF04021U

[Read More](#)

SKU:

Categories:Geriatric Patient Care

Product Description

The KERi™ Basic Mainikin has a non-age-specific appearance.

Complete lifelike range of motion, realistic patient positioning, and non-pinching joints - this Basic Patient simulator moves like a person.

- Lightweight - approximately 28 lbs.
- Full-size adult manikin - measures 58"
- Overall female appearance
- Converts to male with removal of wig and attachment of male genitals

Visual Inspection:

- Normal and Cancerous Mole
- Stage 1 Sacral Ulcer
- Dilated and Constricted pupils

Reddened Skin Folds Patient Care Simulation:

- Bandaging and Wound Dressing
- Bed Baths
- Clothing Changes
- Denture Placement and Removal, Upper and Lower
- Ear Canal Irrigation, Otic Drops, and Hearing Aid Placement
- Eye Irrigation
- Finger and Toe Manipulation with Bandaging - Flexible and Individually Molded
- Hair Care Washing and Combing
- Intramuscular Injection Sites: Arm, Thigh, Buttock
- Oral and Nasal Hygiene: Lavage, Gavage, and Suctioning
- Ostomy Care: Ileostomy and Colostomy Tissue
- Maintenance and Appliance Application, Lavage, and Suctioning
- Patient Positioning
- Patient Transfer Techniques
- Pericare
- Tracheostomy Care, Lavage and Suctioning

	Basic: KERi™ 1020187 GERi™ 1005620	Complete: KERi™ 1013742 GERi™ 1005597	Advanced: KERi™ 1017567 GERi™ 1005606	Auscultation: KERi™ 1020147 GERi™ 1020146
Skin inspection and care	Yes	Yes	Yes	Yes
Bandaging	Yes	Yes	Yes	Yes
Washing	Yes	Yes	Yes	Yes
Denture Placement	Yes	Yes	Yes	Yes
Hearing aid placement	Yes	Yes	Yes	Yes
Ear and eye irrigation	Yes	Yes	Yes	Yes
Ostomy and tracheostomy care	Yes	Yes	Yes	Yes
Positioning and transfer techniques	Yes	Yes	Yes	Yes