



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com


Susie® OB Simulator

Item No. 1013558 [W99999-565]
Weight 11.99 kg
MPN: S500.100.M

[Read More](#)

SKU:

Categories:Obstetrics



Product Description

This versatile skills trainer offers practicing routine delivery techniques and often-seen complications.

Delivered with padded abdominal cover, 3 dilating cervices, 3 vulva inserts, articulating birthing baby, 2 placentas, 3 umbilical cords, 2 umbilical clamps, episiotomy perineum insert, 48 hour postpartum uterine, elevating pillow, silicone lubricant, talcum powder and carrying bag.

Features

- Normal or breech delivery
- Demonstrate internal rotation, expulsion and external rotation
- Inflatable cushion lifts fetus for Leopold maneuver
- Practice intrauterine manipulation to convert breech or transverse to vertex delivery
- Adaptive birth canal to demonstrate dystocia and deal with its relief
- 3 dilating cervices
- 1 fetus for "version" exercises
- 1 padded stomach cover
- 3 Vulva inserts
- 2 placentas
- 3 umbilical cords
- Silicone lubricant and talcum powder
- Instructions
- Carrying bag
- T-shirt and diapers
- Soft carrying bag
- Instruction manual
- Versatile • Portable
- Normal or breech deliveries
- Leopold or breech deliveries
- Intrauterine manipulation
- Cord prolapse
- Placenta previa
- Fundal massage
- Postpartum care

Palpation for Leopold Maneuvers

- Place fetus in vertex, transverse or breech position
- Lower cushion raises fetus for palpation
- Soft stomach cover facilitates palpation or fetal manipulation

Normal Vaginal Delivery

- Birth from ROA or LOP positions
- Demonstrate rotation/ expulsion
- Cervix dilates as labor progresses
- Position placenta to simulate placenta previa
- Practice reduction of umbilical cord

Breech and Vacuum delivery

- Place fetus on birthing mechanism in normal or breech position
- Teach Pinard "leg flip" during breech.
- Pause normal delivery to create need for vacuum augmentation

Postpartum Activity

- Practice fundal massage to firm a "boggy" uterus
- Modules feature medical episiotomy with tears in the labia minora, mediolateral tears, and standard mediolateral episiotomy.

	Birthing Simulator	Original Childbirth Simulator®	3B Birthing Simulator PRO	OB Susie®	Advanced OB Susie®
Item No.	1000002	1005790	1020333	1013558	1019303
Baby					
Fetus	fetus with fontanelles and cranial sutures	male and female articulating fetus	male articulating fetus with fontanelles	articulating fetus	articulating fetus
Palpation of fetal fontanelles	✓	✓	✓	✓	✓
Audible maternal and newborn heart sounds	-	-	-	-	✓
Birth					

Birth mechanism	manual	manual	manual	manual	automatic
Complete, frank, and footling breech births	✓	✓	✓	✓	✓
Cesarean section delivery	-	✓	✓	✓	✓
Complete and incomplete placenta birth	only complete	only complete	✓	only complete	only complete
Artificial Rupture of Membranes (AROM)	-	-	✓	-	-
Intrauterine manipulation	-	✓	✓	✓	✓
Cord prolapse and placenta previa	✓	✓	✓	✓	✓
Maneuvers					
Ritgen's	-	✓	✓	✓	-
Leopold's	✓	✓	✓	✓	✓
Pinard's and Prague's (Breech)	-	✓	✓	✓	✓
Mauriceau-Smellie-Veit (Breech)	-	-	✓	✓	✓
Rubin's and Wood's (Shoulder Dystocia)	-	-	✓	✓	✓
Postpartum					
Episiotomy	-	✓	-	✓	-
Fundal massage	-	✓	-	✓	✓
Boggy uterus		via optional module		✓	✓
Postpartum hemorrhage (PPH)	-	via optional module	via optional module 1021567	-	✓
Clamping and cutting of umbilical cord	-	-	✓	✓	✓