

## **SEM TRAINERS & SYSTEMS**

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile :	+91 88495 63724	Email:	sem@semtrainers.com
Mobile 1:	+91 98791 03905	Website:	www.semtrainers.com



## Deluxe Venous Access Device Model

Item No.	
Weight	
Dimensions	
Brand	
MPN:	
Read More	

1005563 [W43007] 4.05 kg 53 x 38 cm WRS Group 26124

SKU:

Categories: Injections and Punctures







## **Product Description**

Teach patients and healthcare staff the care and use of two different types of catheters with this model. Features left and right subclavian catheter connections that may be flushed with water; implanted venous access system with a "skin" flap over the implant to allow viewing. Made of skinned, lightweight urethane foam, this model may be washed with soap and water, or painted with antiseptics.

Mounted on a collapsible easel and packed in its own carrying case (Replacement door W43008).

<u>Caution!</u> Skin might contain bubbles, please refer to the pictures. However, this does not influence the functionality of the simulator.