



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

I.V. & I.M. Injection Arm

Item No. LF00698U
Weight 9kg
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Injections and Punctures



Product Description

- Realistic Flashback
- Replaceable Skins

- Intramuscular Injection Site complete with Fingerprints

Simple and easy to use, Life/form® Venipuncture and Injection Training Arms are the ideal way to teach venipuncture and injection techniques including starting IVs and introducing Over the Needle IV catheters. This simulator ends your problem of finding willing volunteers, and helps to avoid patient injury and liability while the students are learning new procedures. Realism begins with fine detail - right down to the fingerprints. But most important is the realistic feeling of puncture through the vinyl skin and latex veins. The skin actually rolls as you palpate the vein and the characteristic “pop” can be felt as the needle penetrates the vein. Veins are accessible at the antecubital fossa, along the forearm, and at the back of the hand, making it possible to practice venipuncture at any of the common sites. Practice of intramuscular injection is equally realistic. A foam pad simulates muscle tissue, and the boundaries are identified by natural bonelike landmarks. Life/form® Venipuncture and Injection Arms are designed with replaceable skins and veins. Complete replacement skin kits are available and easy to use. To extend the life of each skin/vein set, an aerosol sealant is available to seal punctures in the tubing and prevent leakage. Includes a Life/form® arm with skin, 22 gauge needle, butterfly set, infusion set, 3 cc syringe, 12 cc syringe, one pint Life/form® blood, fluid supply bags, instruction manual, and hard carry case. Fluid supply stand not included. Size: 35" x 12" x 6". CAUTION: Use only Life/form® Replacements. Do not use Conventional Latex Tubing. This will cause malfunction of the simulator.