



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Geriatric LOR Insert for epidural and spinal injection trainer

Item No. 1020800 [XP61-003]

Weight 0.83 kg

Dimensions 28 x 14 x 7.5 cm

Brand 3B Scientific

[Read More](#)

SKU:

Categories: Injections and Punctures

Product Description

The geriatric LOR (Loss of Resistance) insert for the [Epidural and Spinal Injection Simulator P61](#) lets you practice with the specific anatomical deformations of the spine of an elderly patient.

You can easily replace the standard patient LOR insert (delivered with your Epidural and Spinal Injection Simulator P61) with this optional geriatric insert. It represents the characteristic soft tissue structure (decreased resistance) and the deformed bone structure of a geriatric patient for advanced training possibilities:

Th11-Th12: No access (frustration threshold)

Th12-L1: Centered, normal

L1-L2: Ambilateral, small (partially bony center)

L2-L3: Centered, small (dorsal process L3 shifted to left side)

L3-L4: Centered, partially bony (dorsal process L3 shifted to left side)

L4-L5: Centered, transverse oval

Delivery content:

- 1 Geriatric Insert
- 1 Spinal Tube