



- Fully articulated head, neck, and jaw
- Jaw thrust
- Palpable and visual landmarks
- Palpable brachial pulse (manual)

ECG/Umbilical Cannulation Skin

- Allows practice of umbilical cannulation, cutting, and clamping
- Arrhythmia recognition
- Aspirate fluid through the umbilicus
- Learn proper techniques in umbilical catheterization

Intraosseous Infusion

- Palpable landmarks include the patella, tibia, and tibial tuberosity
- Pressurized system allows aspiration of fluid
- Replaceable bones and skin

IV Arm & Leg

- Realistic flashback from a pressurized system
- Replaceable skin and veins
- Smaller lumen tubing simulates the size of a infant's veins

Manikin can be used with Supraglottic airway devices for airway training.

	Basic Infant CRiSis™	Intermediate Infant CRiSis™	Complete Infant CRiSis™	Deluxe Infant CRiSis™
Item No.	1017248	1021997	1005658	1018146
Airway management	✓	✓	✓	✓
CPR	✓	✓	✓	✓
ECG/Umbilical cannulation Skin	-	-	✓	✓
Intraosseous infusion	-	✓	✓	✓
IV arm	-	✓	✓	✓
IV leg	-	✓	✓	✓
Interactive ECG arrhythmia simulator	optional (1005666)	optional (1005666)	optional (1005666)	included
Carry case	included	included	included	included