



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com

Optional Edema Foot, Geri/Keri

Item No. LF04084U
Weight 0.001kg
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Decubitus Care



Product Description

Add realism while practicing patient assessment skills, including observation of swelling of the foot. Edema simulation can be practiced by pressing on the surface of the skin. The foot uses lifelike simulated tissue pads to demonstrate Stages 1+ through 4+, as well as a non-pitting stage of edema known as brawny edema. Pressing on the surface of the skin may leave an indentation that does not immediately rebound when pressure is released, or it may not indent under pressure, depending on which tissue pad is placed in the foot. A Deep Tissue Injury is also present on the heel. Includes a foot that easily attaches to any GERi™ or KERi™ Patient Care Manikin, five interchangeable pads representing one non-pitting stage and four pitting stages of edema, and an instruction key card.