



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Life/form® GERi™ / KERi™ Manikin IV Arm - Light

Item No. LF04080
Weight
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Injections and Punctures

Product Description

Features replaceable skins and latex veins that actually roll when palpated. A realistic flashback confirms proper needle placement, and the shoulder can be used for intramuscular injection training. For use with all **GERi™** and **KERi™** manikins.