



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Suprapubic Training Model - Light

Item No. SB38897 L

Weight

Make: Nasco Healthcare, USA

[Read More](#)

SKU:SB38897

Categories:Catheterization

Product Description

Become skilled in this increasingly common alternative to urethral indwelling catheters. Feel the normal pressures exerted by the abdominal wall on catheterization. Understand the average length of catheter required to be inserted in order to drain the bladder. Water will flow through catheter when bladder is successfully reached. Lower torso is attached to a vacuum-formed base. Also includes plastic container with lid and fittings, pressure bulb with tubing and connectors, size 16 Foley catheter, instruction manual, and rigid carrying case.