



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Prestan Professional Infant CPR Training Manikin Medium Skin

Item No. PP-IM-100M-MS
Weight
Make: prestan Products, USA

[Read More](#)

SKU:

Categories:BLS Child

Product Description

The PRESTAN Professional Infant Manikin life-like features and weight allows for realistic infant CPR training. Sculpted body contours and realistic skin feel help prepare students for real life rescues. A head that simulates the way an actual infant's head would move. Face/head tilt enhances the CPR learning experience and teaches students how to open the airway during CPR. Equipped with the same CPR feedback features as the Professional Adult and Child Manikins. Fully compliant with current industry guidelines, including the most recent AHA Integrated Feedback Directive. **Prestan Professional Infant CPR Training Manikin (installed with Rate Monitor and Depth Clicker for CPR Feedback and AA Batteries not included). Includes:**









- One Infant Manikin Medium Skin
- 10 Child Face-Shield/Lung-Bags
- Instruction Sheet
- Carry Bag

Warranty Information & Guarantee


PRESTAN offers a THREE YEAR LIMITED WARRANTY against any defects in workmanship or materials under normal use with respect to its PRESTAN TRAINING MANIKINS.

Additional Details:

All PRESTAN Professional Manikins are available with a revolutionary CPR Rate Monitor that allows for instant feedback to both instructors and students. The visual feedback from the CPR Rate Monitor gives students a real life feel of delivering 100-120 compressions per minute. Lights continually signal rate progression, confirming compressions are within the 100- 120 cpm range recommended by the Industry Guidelines. The yellow light will blink if the 120 cpm upper limit is exceeded, warning students to slow down. Instructors can monitor several students quickly and easily. In addition to the visual CPR Rate Monitor, students will also hear a clicker sound as the chest is pushed to the appropriate depth of 2.0 -2.4 inches allowing them to experience the true force needed to deliver real life chest compressions. The CPR Rate Monitor may be purchased separately to upgrade any non-monitored PRESTAN Professional Manikin.

Sleep Mode/No Power	
1st Full Compression (Wake mode/LED Test Mode)	
CPM* < 60	
CPM* 60 – 79	
CPM* 80 – 99	
CPM* 100 - 120	
CPM* > 120	
Compressions Stopped	

*CPM = Compressions per minute

 = Flashing

