



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com


CT QA Phantom JCT II

Item No. PH-54
Weight 3 kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom



Product Description

Renovated CT QA phantom for initial test at the time of reception, as well as follow-up periodical quality check of CT scanners

Features

1. The phantom can be used for initial and follow-up QA tests listed below, described in JIS Z 4752-3-5: 2008 (IEC 61223-3-5: 2004) and Z 4752-2-6: 2012 (IEC 61223-2-6: 2006)
- / 2. Conforming to JIS Z 4923:2015

Training skills / Applications

Set includes

Size (approx.)

Weight (approx.)

Materials

APPLICATIONS

System requirements of software

Update

Axial scan: Slice thickness/spatial resolution/low contrast resolution/noise/mean / Helical scan: Slice thickness HU number/uniformity
1 cylindrical container (fixing screw) / 1 slice thickness unit (axial) / 1 spatial resolution unit / 1 repeated pattern unit / 1 low contrast resolution unit / 2 slice thickness unit (helical) / 1 fixture for the slice thickness unit / 1 fixture for the cylindrical container / 1 Angle adjustment holder (table-top type) / storage case / manual
dia.200×200mm
3kg
Acrylic resin
CT
OS Windows 7 (64 bit), Windows 8.1 Pro (64 bit), memory 4GB, HDD 250 GB
August 10, 2020