



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com


BMD Chart Phantom UHA

Item No. PH-10
Weight 4kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom



Product Description

Bone Mineral Density chart for microdensitometry (MD) method

Features

Training skills / Applications

Set includes

1. 21 steps with different hydroxyapatite content 2. Steps range from 0 to 400 mg/cm², with 20mg/cm² difference each

microdensitometry

1 chart phantom / 1 storage case

Size (approx.)
Packing weight (approx.)
Materials
Update

W3×D21×H1.5cm / W1.2 x D8.3 x H0.6 inch
4kg / 8.8lb
Polyurethane, hydroxyapatite
July 10, 2020