



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Myocardial Phantom HL

Item No. PH-24
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom

Product Description

For the study of high radio accumulation interference in the liver with the myocardial SPECT images

Features

Set includes

Size (approx.)

Packing size (approx.)

Materials

Update

Allows the study of RI liver intake and its effect on the myocardial SPECT Cold defect can be set in the left cardiac muscle Background can be set individually in lung field, mediastinum and right ventricle
1 main Phantom Body / 1 right Lung / 1 left Lung / 1 mediastinum / 1 liver /stomach / 1 heart / 1 work base / 1 screwdriver / 1 petroleum jelly / storage case
W32 x D22 x H31cm / W12.5 x D8.6 x H12.2in
W44 x D39 x H42cm / W17.3 x D15.3 x H16.5in
Main Container: Acrylic Resin / Spine: Epoxy Resin (similar to human in HU) / Heart: Acrylic Resin, Acrylic resin / Lung: Foamed Resin, Water / Screw: Polyacetal Resin
July 7, 2020