



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com




NEW! Ultrasound Guided Tube Feeding Phantom

Item No. US-20
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories:Ultrasound Phantom



Product Description

Anatomically correct neck model for Ultrasound Guided Tube Feeding Training

Features

Training skills / Applications

Set includes

Size (approx.)

Weight (approx.)

Related products

CE certified/CE Certification?

Update

1. Anatomically correct neck model Thyroid Gland Trachea Esophagus Sternum Collarbone Carotid Artery Jugular Vein 2. Realistic neck ultrasound images 3. Usable with real ultrasound diagnostic devices 4. Allows the confirmation of the proper Nasogastric Tube Placement with ultrasound 5. Inner switch for tube insertion direction: trachea/esophagus for tube location image recognition tests (recognition by ultrasound of proper esophagus insertion or wrong tracheal insertion)
Proper handling and manipulation of probe on neck ultrasound Feeding tube insertion Confirmation of tube placement by ultrasound Neck scan to evaluate the tube position.
Main Body Ultrasound neck block Nasogastric Tube Syringe Talcum powder Storage case
W29 x D18 x H17cm(Body) W17 x D10 x H10cm(cube) W11.4 x H6.7 x D7.1inch
1.8kg(body) / 1.1kg(cube)
MW8 Tube Feeding Simulator(NG, OG, and PEG)
No
November 1, 2020