



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com




CVC Insertion Simulator II

Item No. M93UB
Weight 9 kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories:Ultrasound Phantom



Product Description

Conduct a variety of CVC insertion and management simulations with the assistance of 3

interchangeable pads

Features

Training skills / Applications

Set includes

Size (approx.)

Packing size (approx.)

Weight (approx.)

Packing weight (approx.)

Materials

Replacement parts

Optional parts

Related products

Production ? Development Supervision

Notes

Update

1. Realistic ultrasound-guided CVC 2.Landmark puncture pad with anatomically correct vein bifurcations simulates mechanical complications including pneumothorax, mislodging and artery puncture 3.Transparent anatomical block for anatomical understanding and guide wire manipulation

Ultrasound-guided CVC Landmark guided CVC Ultrasound-guided venous access Prevention of mechanical complications

1 male upper torso manikin / 1 landmark puncture pad / 1 ultrasound puncture pad / 1 transparent anatomical block / 1 introductory ultrasound training block / 10 cotton swab for coloring (red) / 10 cotton swab for coloring (blue) / 1 air bulb / 2 plastic jars / 1 irrigation bottle / 1 syringe / 1 sample syringe / 1 sample needle

W40 x D20 x H 34cm / W15.8 x D7.9 x H 13.4in

W52x D45x H38cm/W20x D18x H15in

2.3kg / 5lb

9kg / 19.8lb

Soft resin / Latex free

11347-119 2 landmark puncture pads / 11347-170 2 ultrasound puncture pads / 11347-210 2 introductory ultrasound training blocks / 11347-320 irrigation bottle/11388-400 cotton swab for coloring (red, set of 10)/ 11402-010 cotton swab for coloring(blue, set of 10)/11347-260 1 pulsation unit (for M93UB)

11347-250 pulsation unit for CVC placement pad / 11347-240 2 CVC placement pads

[MW18 Ultrasound-Guided PICC Training Simulator](#) [M93C CVC Insertion Simulator III](#)

Masahiro Tanabe, M.D., Ph.D., Professor Director General Medical Education Center Chiba University School of Medicine / Kinya Sando, M.D., Ph.D., Professor Director Department of Human Dietetics Graduate School of Human Science Osaka Shoin Women's University / Masanori Hoki, M.D., Ph.D., Professor Department of Human Dietetics Graduate School of Human Science Osaka Shoin Women's University / Joho Tokumine, M.D.Ph.D., Department of Anesthesia, Kawatetsu Chiba Hospital

Ultrasound scanner and needles are not included.

October 1, 2020