



**SEM**


## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com




## CVC Insertion Simulator III

**Item No.** M93C  
**Weight** 9 kg  
**Brand** KYOTO KAGAKU, JAPAN

[Read More](#)

**SKU:**

**Categories:**Ultrasound Phantom



## Product Description

CVCIII facilitates a sequence of uninterrupted simulation experience for specialist

## Features

1. Insertion practice is repeatable to obtain the best technique for patient safety.
2. Both landmark guided insertion and ultrasound guided techniques can be practiced.
3. The frictionless tissue-like simulation pad to facilitate the realistic sensation of CVC placement.

Training skills / Applications

Set includes

Size (approx.)

Packing size (approx.)

Weight (approx.)

Packing weight (approx.)

Materials

Replacement parts

Optional parts

Related products

Notes

Update

Perform hand hygiene practices and prepare aseptic field Decide on insertion site via ultrasound guidance or observation and palpation of landmarks Insertion under ultrasound guidance or via landmark technique Confirm flashback Advance guidewire and sheath Insert catheter
1 male upper torso manikin / 2 CVC placement pad / 1 skin for cannulation training / 1 transparent anatomical block / 2 artery tube (transparent): 2 type / 2 vein tube (blue): 2 type / 2 vein pipe / 1 pulsation unit (for landmark technique) / 1 irrigation bottle (for UG technique) / 1 plastic jar / 1 sample needle / each 1 syringe (2.5ml, 50ml) / 1 introductory ultrasound training block REAL VESSELS
W40 x D20 x H34cm / W15.8 x D7.9 x H 13.4in
W52x D45x H38cm/W20x D18x H15in
2.3kg / 5lb
9kg / 19.8lb
Soft resin, hydrogel, hard resin / Latex free
11347-210 2 introductory ultrasound training block REAL VESSELS / 11347-240 2 CVC placement pads / 11347-320 irrigation bottle
11347-260 1 pulsation unit (for M93UB)/11388-400 cotton swab for coloring (red, set of 10)/ 11402-010 cotton swab for coloring(blue, set of 10)/11347-119 2 land mark puncture pads/ 11347-170 2 ultrasound puncture pads
<a href="#">MW18 Ultrasound-Guided PICC Training Simulator</a> <a href="#">M93UB CVC Insertion Simulator II</a>
Ultrasound machine and needles are not included.
October 1, 2020