



**SEM**

## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com






## Nasopharyngeal Swab Collection Simulator

Item No. MW45  
Weight  
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

**SKU:**

**Categories:** Ear, Nose and Throat (ENT)



## Product Description

For training to healthcare professionals who conduct virus testing for the first time

**atures**

1. When swab reaches the right position in the nasal cavity, the swab will be pigmented in blue
2. Ways to use tongue depressor with swabs can be practiced.
3. The shape of nasal cavity between right and left are different.
4. The product comes with cross sectional model

**Training skills / Applications**

**Set includes**

**Size (approx.)**

**Packing size (approx.)**

**Weight (approx.)**

**Packing weight (approx.)**

**Materials**

**Related products**

**[YouTube] JETRO COVID 19 Africa Webinar 2020**

**Update**

Application: Training to healthcare professionals who conduct virus testing for the first time.
1 main model 1 cross sectional model 1 syringe (30ml) 1 disposable cup 10 swabs for coloring 1 storage case 6 replacement pads manual
W42 x D23 x H38 cm W16.8 x D9 x H15.2 in
W49x D29x H53 cm W19x D11.5x H21 in
2.15 kg / 4.73 lb
4kg / 8.8lbs
Soft resin Latex free
<a href="#">Strap-on Intramuscular Injection Trainer</a>
<a href="https://www.youtube.com/watch?v=ljHhgTHRigs">https://www.youtube.com/watch?v=ljHhgTHRigs</a>
December 9, 2024