



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com




Simulator Intravenous ARM III

Item No. MW49
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories:Adult Patient Care



Product Description

Features

Two types of injection pads for venous line placement (forearm/ back of hand) Two types of injection pads for cubital fossa (regular and durable)

Training skills / Applications

Set includes

Size (approx.)

Power source

Power requirements (approx.)

Replacement parts

Related products

Recommended devices

Update

Anatomy of the arm Vacuum tube blood collection Palpation /Intravenous injection IV cannula placement and drip infusion Blood collection with a syringe/Flash-back confirmation
2 arm model / 2 supporting stand / 2 injection pad (R pad) / 2 injection pad (D pad) / 4 injection pad (forearm/antebrachial) / 2 injection pad (back of hand) / 1 circulation pump with a plastic jar / 2 bottle with stand (for sim blood) / 1 simulated blood (10 swab pack) / 2 arm vessels anatomy sleevelet / 1 storage case qty / instruction manual qty
W54 x D11 x H11cm (arm model, each)
AC100V-240V 50/60Hz
5W
11267-300-03 R pad (set of 2) 11430-020 D pad (set of 2) 11388-200 pad for forearm/antebracheal (set of 2) 11388-300 pad for back of hand(set of 2) 11267-300-02 supporting stand 11267-300-05 circualrtion pump 11388-400 simulated blood (10 swab pack) 11430-010 arm vessels anatomy sleevelet
Peripheral Venous Catheter Placement Simulator
IV cannula 23G syringe needle 23G butterfly needle 23G
December 4, 2024