



**SEM**

## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com




## Male Catheterization/Enema Simulator

Item No. MW2A  
Weight 3.5 KG  
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

**SKU:**

**Categories:**Catheterization



### Product Description

High quality simulators for urethral catheterization, perineal care and enema practice

### Features

1. Soft and flexible genitalia allow proper cleaning and procedural preparation
2. Successful insertion confirmation by fluid outflow
3. Enema insertion in lateral position
4. Detachable external genitalia

**Training skills / Applications**

**Set includes**

**Size (approx.)**

**Packing size (approx.)**

**Weight (approx.)**

**Packing weight (approx.)**

**Materials**

**Replacement parts**

**Related products**

**RECOMMENDED DEVICES**

**Update**

|                                                                                                                                                                                                                                                           |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Intermittent catheterization / Indwelling catheterization / Perineal care / Enema / Localization of urethral meatus / Cleansing and aseptic technique / Catheter insertion and placement / Catheter Removal /Manual bladder compression                   |
| 1 lower torso manikin / 1 genital organ unit / 3 valves(one is incorporated in the assembled unit) / 1 irrigation bottle / 1 rectum tube / 3 irrigation drainage tubes / 1 lubricant /1 storage case / 1 instruction manual                               |
| W34 x D33 x H24 cm / W13.5 x D13 x H9.5 in                                                                                                                                                                                                                |
| W52 x D45 x H31cm / W20.5 x D17.7 x H12.3in                                                                                                                                                                                                               |
| 3.5kg / 7.7Ibs                                                                                                                                                                                                                                            |
| 14.5kg / 32Ibs                                                                                                                                                                                                                                            |
| Soft resin / Latex free                                                                                                                                                                                                                                   |
| 11381-100-01 clip (small) 11381-100-02 clip (large) 11381-100-03 urinary bladder reservoir 11381-100-04 rectum unit 11381-300 6 replacement valves 11381-800 lubricant 11381-410 replacement male skin 11381-510 replacement male meatus for MW2A(a pair) |
| <a href="#">MW52 Female Catheterization/Enema Simulator</a><br><a href="#">MW53/54/55 Strap-on Catheterization/Enema Simulator</a><br><a href="#">MW58/59 Suprapubic Catheter Care Simulator</a> <a href="#">MW23 Adult diaper training simulator</a>     |
| Urinary catheter: 16Fr                                                                                                                                                                                                                                    |
| October 20, 2020                                                                                                                                                                                                                                          |