



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com


Peripheral Venous Catheter Placement Simulator

Item No. MW9
Weight 0.9kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Intravenous (I.v.) and Arterial



Product Description

Holistic simulation of peripheral intravenous access management in just one arm

Features

1. IVC insertion and management from start to finish
2. Two puncture sites: median antebrachial vein and dorsal veins
3. Durable puncture pad

Training skills / Applications

Set includes

Size (approx.)

Packing size (approx.)

Weight (approx.)

Packing weight (approx.)

Replacement parts

Related products

Recommended Devices

Update

Use of tourniquet / Identifying puncture site / Sterilization of puncture site / Puncture with IV cannula / Confirmation of flashback / Application of pressure on vein / Preparing and administering IV medications (simulated) / Dressing and securing cannula
1 arm model / 2 injection pads for dorsal veins of hand / 2 injection pads for median antebrachial vein / 1 bottle for simulated blood / 1 stand for bottle of simulated blood / 1 infusion set / 1 infusion bag / 1 syringe (50ml) / 1 plastic beaker / 10 simulated blood sticks (swab type) / 1 carrying bag / 1 instruction manual
W52 x D10 x H10cm / W20.5 x D3.9 x H3.9in
W68x D23x H20cm/W27x D9x H8in
0.9kg / 2lb
2.6kg/5.7lb
11388-200 2 injection pads for median antebrachial vein / 11388-300 2 injection pads for dorsal veins of hand / 11388-400 10 simulated blood sticks (swab type)
MW50 Strap-on Venipuncture Trainer
IV catheter needle 22G, butterfly needle 23G
December 4, 2024