



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com


Ultrasound-Guided Pericardiocentesis Simulator

Item No. MW15
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories:Ultrasound Phantom



Product Description

Practice patient safety during pericardiocentesis with ultrasound guidance / This simulator allows trainees to insert the needle under ultrasound guidance, pierce the "pericardial sac" and aspirate pericardial fluid.

Features

Training skills / Applications

Set includes

Size (approx.)

Packing size (approx.)

Weight (approx.)

Packing weight (approx.)

Materials

Replacement parts

Optional parts

Related products

Production ? Development Supervision

Recommended Devices

Notes

Update

1. Durable and replaceable puncture pad 2. Practice both the subxiphoid approach and parasternal approach through palpable and ultrasonic landmarks. 3. Realistic needle tip sensation during puncture of the "pericardial sac"

Patient positioning / Visualization of pericardial fluid under ultrasound scanning / Landmark palpation / Needle insertion to the pericardial space / Aspiration of pericardial effusion

1 adult chest model / 1 pericardiocentesis unit / 1 pillow for positioning / 1 instruction manual / 1 syringes (50ml) / 1 joint hose

Torso size: W38 x D25 x H48 cm , 15 x 9.9 x 18.9 in / Pericardio pad size: W16 x D14 x H21 cm , 6.3 x 5.6 x 8.3 in

W65 x D48 x H48 cm / W25.5 x D19 x H19 in

9.65kg / 4.4lb

12kg / 26.5lb

Soft resin / Latex free

11394-010 pad for MW15

11383-010 mid axially line unit puncture pad / 11383-020 mid-scapular line unit puncture pads / 11383-030 replacement lung

[MW17 Ultrasound-Guided Thoracentesis/ Pericardiocentesis Simulator](#) [MW4 Ultrasound-Guided Thoracentesis Simulator](#) [M4A Ultrasound-Guided Thoracentesis Simulator - Strap-on set](#) - [MW15 Ultrasound-Guided Pericardiocentesis Simulator](#)

Takahiro Amano, M.D., Vice President, Senior Vice-Dean, Postgraduate School, Professor and Director, Center of Postgraduate Medical Education, International University of Health and Welfare Honorary Director, Sanno Medical Center

Pericardiocentesis: 18 G needle

Needles are not included.

October 20, 2020