



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Gluteus Maximus Injection Simulator

Item No. M74
Weight 3.2kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Injections and Punctures



Product Description

This simulator is designed to practice in the three ways of injection into the buttocks.

Features

Training skills / Applications

Set includes

Size (approx.)

Weight (approx.)

Power source

Power requirements (approx.)

Replacement parts

CE certified/CE Certification?

Update

1. The model is made from three individual layers, which are skin, muscular and skeleton. Therefore it allows identifying the injection spot as well as practicing with a patient. 2. The model is designed for three methods: Quarter-Trichotomy, Hochstetter Method and Clark Method. 3. Attached indicator shows the accuracy of injection. If the injection spot is correct, the syringe can inject liquid.

Injection into the buttocks (Quarter-Trichotomy, Hochstetter Method and Clark Method)

1 main body / 1 stand for Hochstetter method / 1 controller / 1 set of sensor pads / instruction manual

W32 x D35 x H20cm

3.2kg

(AC100V 50/60Hz)

2W

11213-010 Skin / 11213-020 Sensor for Quarter-Trichotomy / 11213-030 Sensor for Hochstetter Method / 11213-040 Sensor for Clark Method / 11213-050 Controller

No

October 10, 2020