



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Female Condom Model

Item No. 1005826 [W45154]

Weight 1.64kg

Dimensions 17.8 x 24 x 14 cm

MPN: S509

[Read More](#)

SKU:

Categories:Condom Trainers

Product Description

This model represents an anteverted uterus in a simulated pelvic cavity with soft vulva and vagina and a soft plastic stomach cover. A perfect tool to demonstrate the use of a female condom, contraceptive sponge and cervical cap. Supplied with carrying bag.

Female Condom Model (W45154) — Anteverted Uterus Pelvic Model for Female Barrier Method Training

The Female Condom Model (Item No. 1005826 [W45154], 1.64 kg) represents an anteverted uterus in a simulated pelvic cavity, complete with soft vulva and vagina and a soft plastic stomach cover — giving students a realistic, anatomically accurate platform for female barrier contraception training.

A versatile training tool for multiple barrier methods. This model is built to demonstrate the correct use of a female condom, contraceptive sponge, and cervical cap — covering three distinct female barrier contraception methods with a single teaching model, rather than requiring separate models for each.

Soft, realistic materials for genuine hands-on practice. The soft vulva, vagina, and stomach cover give students a tactile, realistic experience when practicing insertion technique — important for building genuine confidence before working with real patients.

Comes ready for classroom or field use, supplied with a carrying bag for easy transport and storage between training sessions.

Why Buy from SEM Trainers? SEM Trainers & Systems has over 30 years of experience supplying medical education resources to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What contraceptive methods can be practiced on this model? A: Female condom, contraceptive sponge, and cervical cap insertion technique.

Q: What's included? A: The pelvic model itself and a carrying bag for transport and storage.

Q: Is this suitable for both classroom and clinical training settings? A: Yes — its realistic materials and portable carrying bag make it suited to both academic and field-based reproductive health training.

Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.