



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Standard Breast Self Examination Model, Beige

Item No. 1017889 [W43004]
Weight 0.9kg
Dimensions 25 x 7.5 x 15 cm
Brand WRS Group
MPN: 26500

[Read More](#)

SKU:

Categories: Women's Health Education



Product Description

This soft BIOLIKE 2™ model mimics the feel of real tissue and features palpable lumps that teach fingers how to locate potential abnormalities. The model's non-palpable lumps can be seen only with the included penlight, underscoring the importance of mammography. A great tool for instruction in breast self-exam and breast self-awareness programs. **Standard Breast Self Examination Model, Beige (W43004)**

The Standard Breast Self Examination Model, Beige (Item No. 1017889 [W43004], 0.9 kg) by WRS Group is a soft BIOLIKE 2™ model that mimics the feel of real tissue, featuring palpable lumps that teach fingers how to locate potential abnormalities.

Two lump types teach two important lessons. Palpable lumps train hands-on detection technique directly, while the model's non-palpable lumps can only be seen with the included penlight — underscoring the importance of mammography even when a self-exam feels clear.

A great tool for instruction in breast self-exam and breast self-awareness programs, giving both patients and trainees hands-on practice building the tactile skill self-examination requires.

Why Buy from SEM Trainers? SEM Trainers & Systems has over 30 years of experience supplying medical education resources to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Frequently Asked Questions

Q: What's the difference between the palpable and non-palpable lumps? A: Palpable lumps can be felt by hand to train detection technique; non-palpable lumps can only be seen with the included penlight, reinforcing why mammography matters even after a normal self-exam.

Q: What material is the model made from? A: Soft BIOLIKE 2™ material designed to mimic the feel of real breast tissue.

Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.