



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Suture Practice Leg, Light

Item No. LF01034U
Weight 1.087kg
Dimensions 62 x 18 x 13 cm
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Suturing and Bandaging

Product Description

Made with a soft vinyl skin over a core of stitchable foam to provide a lifelike suturing experience for students or a realistic suturing demonstration by an instructor. Realistic skin texture with wrinkles and pores. Soft and pliable for easy sewing, the skin is tough enough that sutures will not pull out when tightened. Each unit provides three "wounds." These "wounds" can be sutured repeatedly until the skin around them is finally worn out; then, new "wounds" can be cut. Although disposable, the trainers provides hundreds of suturing experiences before wearing out. Conservatively, over one hundred cuts can be made on the trainer and each of these cuts can be sutured several times. Suture Practice Leg is

complete with a starter suturing kit and instruction booklet.

Why Choose the Life/form® Suture Practice Leg, Light (LF01034U) for Your Simulation Lab in India?

The Suture Practice Leg, Light (Item No. LF01034U) by Nasco Healthcare, USA is a lower-limb suture trainer constructed from soft vinyl skin over a stitchable foam core — a material combination that provides realistic needle-passage resistance and natural suture tension for both student practice and instructor demonstration. The leg is available in a light skin tone. For MBBS surgery skills labs, nursing colleges, paramedic training programmes, and hospital in-service courses in India that require a dedicated lower-limb suture trainer, this leg provides the anatomical form and skin fidelity of the leg that arm-based trainers cannot replicate — important for training on wound types commonly encountered in road traffic injuries, falls, and lower-limb surgical procedures.

The skin surface features realistic wrinkles and pores for a lifelike appearance, and is soft and pliable for easy needle passage while remaining tough enough that sutures do not pull out when tightened — a common failure mode in lower-quality suture pads that creates poor tension habits in trainees. The leg comes with three pre-formed wounds ready for immediate use. Once the skin around a wound site is finally worn, new wounds can be cut in fresh locations; conservatively, over one hundred cuts can be made across the trainer and each cut can be sutured several times, providing hundreds of suturing experiences in total from a single unit. A starter suturing kit and instruction booklet are included.

Why Buy from SEM Trainers?

SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support.

Frequently Asked Questions

Q: What does "Light" in the product name refer to? A: It refers to the skin tone of the trainer — this is the light complexion variant of the Suture Practice Leg.

Q: How many suturing sessions does this leg provide? A: The leg comes with three pre-formed wounds. Conservatively, over 100 cuts can be made across the trainer and each cut can be sutured several times, providing hundreds of suturing experiences from a single unit.

Q: Does a suturing kit come included? A: Yes. A starter suturing kit and instruction booklet are included.

Q: Why train on a leg rather than a flat suture pad? A: The cylindrical three-dimensional anatomy of the leg more closely replicates the tissue geometry and needle angle required for lower-limb wound closure — particularly relevant for trauma and emergency surgery training.