



## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724

**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com

**Website:** www.semtrainers.com



## SAM Student Auscultation Manikin Storage/Carry Case w/ wheels and handle

Item No. 1020109  
Weight 0.001kg  
Brand Cardionics  
MPN: 718-8805

[Read More](#)

**SKU:**

**Categories:**Auscultation

## Product Description

### SAM Student Auscultation Manikin Storage/Carry Case — Wheeled Transport for Simulation Labs in India

The SAM Student Auscultation Manikin Storage/Carry Case (Item No. 1020109, MPN: 718-8805) by Cardionics, USA is a wheeled storage and transport case built for the SAM family of student auscultation manikins — SAM 3G, SAM II, SAM Basic, and PAT. With integrated wheels and a pull handle, it's designed for institutions that move manikins between skills labs, exam halls, or storage rooms regularly, without needing two people to lift and carry the unit each time.

For nursing colleges, medical colleges, and AYUSH institutions running auscultation training across multiple classrooms or shared skills labs, a dedicated case protects the manikin's internal electronics and sound system during transport and reduces handling damage between sessions — a common wear point when trainers are shared across departments or moved for exams and demonstrations.

**Why Buy from SEM Trainers?** SEM Trainers & Systems is an authorised distributor of Cardionics products in India with over 30 years of experience supplying medical simulators and auscultation training accessories to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Q: Which SAM manikins does this case fit? A: It's built for the SAM family of student auscultation manikins (SAM 3G, SAM II, SAM Basic, PAT). Confirm your specific manikin model with our team before ordering to ensure fit.

Q: Why use a wheeled case instead of carrying the manikin directly? A: It protects the manikin's internal speaker and electronics during transport and makes it easier for one person to move the unit between rooms — useful for institutions sharing a manikin across multiple labs or classes.

Contact us at [sem@semtrainers.com](mailto:sem@semtrainers.com) or call +91-88495 63724 for pricing, dimensions, and availability.