



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Basic Murmurs Sound Menu

Item No. 1018176 [W49424]
Weight 0.07kg
Dimensions 16 x 7 x 2 cm
MPN: ST2M-142 - Murmurs

[Read More](#)

SKU:

Categories:Auscultation

Product Description

Clean, genuine content this time — no duplication, no obvious errors. Companion card to the two you've already got covered for the same base unit.

[Cert] Item No. 1018176 [W49424], MPN ST2M-142 – Murmurs, 0.07 kg, 16 × 7 × 2 cm, compatible with **W49407** (same unit as the ACLS, Aortic Regurgitation, and W49423 Basic Heart Sounds cards). No placeholder image issue flagged this time beyond the missing product photo itself, and the sound list is internally coherent: it's organized by murmur *timing* (holosystolic, mid-systolic, late systolic, early diastolic, mid-diastolic, holodiastolic, combined presystolic patterns, continuous) rather than by named pathology — a genuinely different teaching approach from the Aortic Regurgitation card, which was organized by diagnosis.

That timing-based structure is worth foregrounding in the copy, since it's the actual pedagogical differentiator versus your other cardiac sound menu cards.

DESCRIPTION (copy-paste into WordPress):

Why Choose the Basic Murmurs Sound Menu for Your Simulation Lab in India?

The Basic Murmurs Sound Menu (Item No. 1018176 / W49424, MPN: ST2M-142, 0.07 kg, 16 × 7 × 2 cm) by Nasco Healthcare, USA is a 16-sound auscultation card for the W49407 simulator, built around murmur **timing** rather than named diagnoses — the core skill students need before they can correctly localize a murmur to a specific valve lesion.

The card progresses through holosystolic murmurs (including crescendo and late-crescendo variants), mid- and late-systolic murmurs, early diastolic murmurs (long blowing, brief, and musical variants), mid-diastolic murmurs with and without presystolic accentuation, holodiastolic murmurs, and a continuous murmur — giving students a structured framework for classifying any murmur by its position in the cardiac cycle before attempting to name the underlying lesion.

This timing-first approach complements diagnosis-focused cards like the Aortic Regurgitation Sounds Menu: students who master timing classification here are better equipped to work through named-pathology cards with real clinical reasoning rather than pattern memorization.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and auscultation training accessories to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Q: How is this different from a diagnosis-specific card like Aortic Regurgitation? A: This card organizes sounds by *when* in the cardiac cycle the murmur occurs (systolic, diastolic, continuous), building the timing-classification skill needed before tackling diagnosis-specific cards.

Q: Which simulator is this compatible with? A: The W49407 auscultation simulator by Nasco Healthcare. Confirm compatibility before ordering.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.

Basic Murmurs Sound Menu - Auscultation Sounds for W49407

1. Holosystolic - blowing crescendo
2. Holosystolic
3. Holosystolic - late crescendo
4. Mid systolic
5. Late systolic
6. Late systolic
7. Late systolic - whoop
8. Early diastolic - long blowing
9. Early diastolic - brief
10. Early diastolic - musical
11. Mid diastolic & rumble
12. Holodiastolic
13. Mid diastolic & presystolic
14. Long mid systolic & presystolic
15. Holosystolic & mid diastolic
16. Continuous