



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com

Paced Beats Sound Menu

Item No. 1018189 [W49416]

Weight 0.07kg

Dimensions 16 x 7 x 2 cm

MPN: ST2M-49

[Read More](#)

SKU:

Categories: Auscultation



Product Description

Why Choose the Paced Beats Sound Menu for Your Simulation Lab in India?

The Paced Beats Sound Menu (Item No. 1018189 / W49416, MPN: ST2M-49, 0.07 kg, 16 × 7 × 2 cm) by Nasco Healthcare, USA is a specialized auscultation menu card for the W49407 auscultation simulator, focused entirely on recognizing pacemaker function and malfunction patterns — a competency general cardiology sound menus don't cover and one that's increasingly relevant as pacemaker implantation rates rise across Indian cardiac care centres.

The 16 tracks are organized by pacing mode: **ventricular pacing** (tracks 1–4: full capture, intermittent capture, and pacing during atrial fibrillation and atrial flutter), **demand pacing** (tracks 5–9: atrial paced,

PVCs, and pacing complicated by atrial fibrillation, right and left bundle branch block), and **AV sequential pacing** (tracks 10–15: PVCs, full capture, occasional capture, and intermittent capture) — with a final ventricular intermittent-capture track (16) for direct comparison against the earlier ventricular series. This mode-by-mode structure trains students to recognize not just that a pacemaker is present, but *which mode* it's operating in and whether it's capturing reliably — a distinction with direct clinical consequences for device follow-up.

For cardiology nursing programmes, ICU/CCU training centres, and postgraduate cardiac care courses in India where pacemaker patient assessment is part of the curriculum, this card fills a gap most general auscultation training sets don't address.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and auscultation training accessories to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Q: What pacing modes does this card cover? A: Ventricular pacing, demand pacing, and AV sequential pacing — including both reliable capture and capture failure/intermittent capture patterns for each mode.

Q: Who should train on this card? A: Cardiology nursing staff, ICU/CCU teams, and postgraduate cardiac care trainees who need to recognize pacemaker function on physical exam, not just on ECG.

Q: Which simulator is this compatible with? A: The W49407 auscultation simulator.

Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.

Paced Beats Sound Menu - Auscultation Sounds for W49407 1. Ventricular, 100% capture 2.

Ventricular, intermittent 3. Ventricular, atrial fibrillation 4. Ventricular, 100% atrial flutter 5. Demand,

atrial paced 6. Demand PVCs 7. Demand, atrial fibrillation, PVCs, RBBB 8. Demand, intermittent,

LBBB 9. Demand, atrial fibrillation 10. AV Sequence, PVCs 11. AV Sequence, 100% capture 12. AV

Sequence, 100% capture 13. AV Sequence, 100% capture 14. AV Sequence, occasional capture 15.

AV Sequence, intermittent capture 16. Ventricular, intermittent capture