



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Advanced Airway Larry Trainer Head, Stand

Item No. LF03685U

Weight 12kg

Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Airway Management Adult

Product Description

Life/form® Advanced “Airway Larry” Airway Management Trainer offers tongue swelling and laryngospasm in addition to all the features on the standard models. Your students can now be presented with the additional challenges they may face in the real world. This Life/form® airway management trainer simulates a nonanesthetized patient for practicing intubation, ventilation, and suction techniques. Realistic anatomy and landmarks include teeth, tongue, oral and nasal pharynx, larynx, epiglottis, arytenoids, false cords, true vocal cords, trachea, lungs, esophagus, and stomach. The trainer allows you to practice oral, digital, and nasal intubation, as well as E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and KING System insertion. Applying pressure to the cricoid cartilage changes the position of the trachea and closes the esophagus, allowing realistic practice of the Sellick maneuver. With its slightly anterior position, swelling tongue, and vocal cords, the Advanced “Airway Larry” is a great trainer for introductory as well as advanced training. Suction techniques and proper

cuff inflation may also be practiced and evaluated. Features durable, rugged, one-piece construction and bifurcated lungs. Head is mounted on a sturdy stand. Includes pump spray lubricant and hard carrying case.

The Life/form® Advanced Airway Larry Trainer Head with Stand (Item No. LF03685U / 1019249) by Nasco Healthcare, USA is a 12 kg advanced airway management training head mounted on a sturdy stand — the stand-format configuration in the Advanced Airway Larry range, with a hard carry case included. From the family comparison table, it is the only Advanced Airway Larry head configuration that combines laryngospasm, tongue oedema, digital intubation, suctioning, all advanced airway devices, and manual carotid pulse simultaneously. The comparable board-mounted Advanced Trainer Head (1020961) has laryngospasm, tongue oedema, stomach rise, and all advanced airway devices, but lacks the carotid pulse. This configuration has the carotid pulse but no stomach rise — making it the correct choice for ACLS and advanced airway programmes where integrated pulse assessment is a training requirement alongside difficult airway management.

The slightly anterior airway position combined with laryngospasm and tongue oedema simulation recreates the genuinely challenging airway geometry that trainees encounter in uncooperative, sedated, or critically ill patients — a training condition absent from the standard Airway Larry head (1012740). Digital, oral, and nasal intubation are all supported, alongside the full advanced airway device set: ET, EOA, PTL, LMA, EGTA, Combitube®, and King Systems. The Sellick manoeuvre is practised by applying pressure to the cricoid cartilage, which changes the tracheal position and closes the oesophagus — providing correct tactile feedback for this critical airway protection technique. Suction technique and correct cuff inflation can both be practised and evaluated. The durable one-piece construction with bifurcated lungs provides reliable bilateral auscultation confirmation of correct tube placement. The stand mounting presents the head at the correct approach angle for standing intubation technique training, distinct from the supine board-mounted configurations. A manual is available for download on the product page.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support. All products carry genuine manufacturer warranty and meet ISO 9001 quality standards.

Frequently Asked Questions

Q: How does this configuration (LF03685U) differ from the other Advanced Airway Larry head (1020961)? A: LF03685U is stand-mounted and includes a manual carotid pulse but no stomach rise. 1020961 is board-mounted, includes stomach rise, but does not have a manual carotid pulse. Both have laryngospasm, tongue oedema, and the full advanced airway device set. Choose LF03685U when integrated pulse assessment training is a programme requirement.

Q: Is the hard carry case included? A: Yes. A hard carry case and pump spray lubricant are included as standard with LF03685U.

Q: What is the warranty and how do I get support? A: All Nasco Healthcare products carry standard manufacturer warranty. For support, reach us at sem@semtrainers.com or call +91-88495 63724.

	Head/torso	Ideal for	Digital intubation	Laryngospasm & Tongue edema	Stomach rise	Sucti-oning	Jaw thrust	Manual carotid pulse	E.T., E.O.A., P.T.L, L.M.A., Combitube®, and King System insertion	Additional features
Adult Airway Management Trainer Head Item No: 1012740	male head	introductory classes	Yes	-	-	Yes	-	Yes	Yes	can be attached to any of the adult CRiSis® Manikins
Advanced Trainer Head Item No: 1020961	male head mounted on a board	introductory and advanced classes	Yes	Yes	Yes	Yes	-	-	Yes	can be attached to any of the adult CRiSis® Manikins
Adult Airway Management Trainer Item No: 1005633	male head mounted on a board	testing purposes and Combitube® practice	Yes	-	-	Yes	Yes	-	Yes	can be attached to any of the adult CRiSis® Manikins
Economy Adult Airway Management Trainer Item No: 1005780	male head mounted on a board	rescue & advanced airway management training	-	-	-	Yes	-	-	only LMAs and Combitube®	durable construction and bifurcated lung

Adult Airway Management Trainer <i>Item No:</i> 1017871	torso with a female head	difficult airway and CPR training	Yes	-	-	Yes	Yes	Yes	Yes	durable construction and bifurcated lung
Advanced Trainer Head <i>Item No:</i> 1019249	male head mounted on a board	advanced airway management training	Yes	Yes	-	Yes	-	Yes	Yes	durable construction and bifurcated lung
S.A.L.A.D. Simulator <i>Item No:</i> 1021583	male head	suction assisted laryngoscopy and airway decontamination	Yes	-	-	Yes	-	Yes	Yes	rigid suction catheter distracts lower mandible and tongue to permit laryngoscope insertion
Airway Management Trainer <i>Item No:</i> 1019928	male torso	airway management and CPR training	Yes	-	-	Yes	-	Yes	Yes	durable construction and bifurcated lung
Adult Deluxe Airway Management Trainer <i>Item No:</i> 1018851	male head mounted on a board	cricothyrotomy and advanced airway management	-	Yes	Yes	-	-	-	-	lifelike skin that will not be damaged by aggressive intubation
Adult Airway Management Trainer <i>Item No:</i> 1005774	male torso mounted on a board	rescue airway and anatomy training	-	Yes	Yes	-	-	-	-	review of the pulmonary system and oral/nasal pharyngeal spaces