



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Basic GERi™ Nursing Skills Elderly-Care Manikin

Item No. 1005620 [W44077]
Weight 10.785kg
Dimensions 94 x 52 x 28 cm
MPN: LF04040U

[Read More](#)

SKU:

Categories:Male Examination



Product Description

Why Choose the Basic GERi™ Nursing Skills Elderly-Care Manikin for Your Nursing College or Skills Lab in India?

The Basic GERi™ Nursing Skills Elderly-Care Manikin (Item No. 1005620 / W44077, MPN: LF04040U) is a full-size, lifelike elderly-care simulation manikin from Nasco / Lifecast Body Simulation, USA — purpose-built for nursing colleges, ANM/GNM programmes, BSc Nursing and post-basic nursing institutes, and hospital in-service training departments in India where geriatric patient care competency is a core curriculum requirement. At 10.785 kg and 147 cm (58 inches), the GERi™ is a realistic, manageable training body that moves like a person — with a full, non-pinching range of motion that allows trainees to practice safe patient handling exactly as they would on a ward.

The GERi™ features an authentic elderly appearance with skin wrinkles and folds, making it the only nursing manikin that prepares students for the specific visual cues and handling considerations of older patient care. The manikin defaults to female presentation and converts to male with removal of the wig and attachment of the included male genitalia — making it suitable for both female and male geriatric patient care training without the need for two separate manikins.

Visual Inspection Skills: Four distinct visual findings are presented on the body for skin assessment training — a normal mole and a cancerous mole (for lesion differentiation), a Stage 1 sacral pressure ulcer (for early-stage pressure injury recognition), dilated and constricted pupils (for neurological and medication-effect assessment), and reddened skin folds (for moisture-associated skin damage identification).

Patient Care Simulation — What the Basic GERi™ Trains:

The Basic tier covers the foundational nursing skills module required in all Indian nursing curricula:

- **Bandaging and Wound Dressing**
- **Bed Baths** — full body wash simulation
- **Clothing Changes** — with realistic patient resistance and positioning
- **Denture Placement and Removal** — upper and lower sets
- **Ear Canal Irrigation, Otic Drops, and Hearing Aid Placement**
- **Eye Irrigation**
- **Finger and Toe Manipulation with Bandaging** — individually moulded, flexible digits
- **Hair Care** — washing and combing
- **Intramuscular Injection** — arm, thigh, and buttock sites
- **Oral and Nasal Hygiene** — lavage, gavage, and suctioning
- **Ostomy Care** — ileostomy and colostomy tissue, appliance application, lavage, and suctioning
- **Patient Positioning**

- **Patient Transfer Techniques**
- **Pericare**
- **Tracheostomy Care** — lavage and suctioning

What the Basic GERi™ Does Not Include (available in higher tiers): Gastrostomy, enema, pap smear, urinary catheterisation, nasogastric tube placement, IV insertion, blood pressure monitoring, and heart/lung auscultation are available in the Complete, Advanced, and Auscultation versions of the GERi™ — all available from SEM Trainers.

Ideal for OBRA-aligned training programmes and OSCE nursing station setups. Note: carry case is not included.

Why Buy from SEM Trainers? SEM Trainers & Systems has over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states in India. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support

Q: Is this a male or female manikin? A: The Basic GERi™ ships with a female presentation by default. It converts to male by removing the wig and attaching the included male genitalia — no additional purchase required.

Q: What is the height and weight of the manikin? A: The GERi™ measures 147 cm (58 inches) tall and weighs 10.785 kg. Dimensions of the shipping package are 94 × 52 × 28 cm.

Q: Does this manikin support catheterisation or NG tube placement? A: No — those skills are available in the Complete GERi™ (1005597). The Basic tier covers the core nursing skills module: bandaging, bed baths, injections, ostomy, tracheostomy, pericare, positioning, and transfer. See the comparison table below for the full feature breakdown across all four GERi™ tiers.

Q: Can the GERi™ be used for OSCE station rotation? A: Yes. At 10.785 kg with realistic range of motion and quick convertibility between male and female presentation, the GERi™ is well suited to OSCE nursing skills stations, in-service ward training, and faculty demonstration.

Q: Is a carry case included? A: No. The carry case is not included with the Basic GERi™. It is available separately — contact SEM Trainers for availability.

Q: What is the difference between the Basic, Complete, Advanced, and Auscultation versions?

A: All four versions share the core nursing skills module (skin inspection, bandaging, bed baths, injections, ostomy, tracheostomy, pericare, positioning, and transfer). The Complete adds catheterisation, NG tube, enema, gastrostomy, and pap smear. The Advanced adds IV insertion and blood pressure monitoring. The Auscultation version adds heart and lung sound recognition. All four are available from SEM Trainers.

GERi™ has an elderly appearance with skin wrinkles and folds. Complete lifelike range of motion, realistic patient positioning, and non-pinching joints - moves like a person. • Full-size adult manikin - measures 58". • Lightweight - approximately 28 lbs. • Overall female appearance • Converts to male with removal of wig and attachment of male genitals Visual Inspection: • Normal and Cancerous Mole • Stage 1 Sacral Ulcer • Dilated and Constricted pupils • Reddened Skin Folds Patient Care Simulation: • Bandaging and Wound Dressing • Bed Baths • Clothing Changes • Denture Placement and Removal, Upper and Lower • Ear Canal Irrigation, Otic Drops, and Hearing Aid Placement • Eye Irrigation • Finger and Toe Manipulation with Bandaging - Flexible and Individually Molded • Hair Care Washing and Combing • Intramuscular Injection Sites: Arm, Thigh, Buttock • Oral and Nasal Hygiene: Lavage, Gavage, and Suctioning • Ostomy Care: Ileostomy and Colostomy Tissue • Maintenance and Appliance Application, Lavage, and Suctioning • Patient Positioning • Patient Transfer Techniques • Pericare • Tracheostomy Care, Lavage and Suctioning Ideal for OBRA Required Training. Carry case not included.

	Basic: KERi™ 1020187 GERi™ 1005620	Complete: KERi™ 1013742 GERi™ 1005597	Advanced: KERi™ 1017567 GERi™ 1005606	Auscultation: KERi™ 1020147 GERi™ 1020146
Skin inspection and care	Yes	Yes	Yes	Yes
Bandaging	Yes	Yes	Yes	Yes
Washing	Yes	Yes	Yes	Yes
Denture Placement	Yes	Yes	Yes	Yes
Hearing aid placement	Yes	Yes	Yes	Yes
Ear and eye irrigation	Yes	Yes	Yes	Yes
Ostomy and tracheostomy care	Yes	Yes	Yes	Yes
Positioning and transfer techniques	Yes	Yes	Yes	Yes
Pericare	Yes	Yes	Yes	Yes
Intramuscular injections	Yes	Yes	Yes	Yes
Gastrostomy	-	Yes	Yes	Yes
Pap smears and douching	-	Yes	Yes	Yes
Enema administration	-	Yes	Yes	Yes
Prostate exam	-	Yes	Yes	Yes
Urinary catheterisation	-	Yes	Yes	Yes

Nasogastric tube placement	-	Yes	Yes	Yes
IV insertion	-	-	Yes	
Blood pressure monitoring	-	-	Yes	
Recognition of heart sounds	-	-	-	Yes
Recognition of lung sounds	-	-	-	Yes