



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



CPR Monitor

Item No. 1005800 [W45050]

Weight 0.723kg

MPN: S180

[Read More](#)

SKU:

Categories:BLS and CPR Accessories

Product Description

for W45117 Displays cadence and depth of cardiac compression and airway ventilation.

Brand identification needed. Item No. 1005800 [W45050], 0.723 kg, MPN S180. Make field absent. The S-prefix MPN (S180) and W45050 catalogue number point to **Gaumard Scientific, USA** — the S-series MPN is Gaumard's nomenclature, consistent with all other Gaumard products in this session. Compatible with W45117 — [Likely] a Gaumard CPR manikin. Displays cadence and depth of cardiac compression and airway ventilation — this is a standalone CPR feedback monitor unit, not integrated into a manikin. Product image is WooCommerce placeholder.

DESCRIPTION (copy-paste into WordPress product description):

Why Choose the Gaumard CPR Monitor for Your Simulation Lab in India?

The Gaumard CPR Monitor (Item No. 1005800 / W45050, MPN: S180) by Gaumard Scientific, USA is a 0.723 kg standalone CPR feedback monitoring unit compatible with the Gaumard CPR manikin W45117 — designed to display real-time cadence and depth of cardiac compression and airway ventilation during CPR training sessions. For BLS certification centres, PALS training centres, nursing institutes, and hospital in-service programmes in India running CPR training on the compatible W45117 manikin, this monitor provides the instructor-visible real-time quality feedback that post-session verbal debrief cannot replace.

The CPR Monitor displays the two most critical CPR quality parameters simultaneously: compression cadence (rate, in compressions per minute) and compression depth — the two variables that AHA guidelines specify as the primary determinants of CPR quality. When both are visible to the instructor during a session, incorrect technique can be corrected in real time during the compression cycle rather than after the trainee has stopped, which is the learning moment where correction has the greatest impact on technique calibration. Airway ventilation is also monitored and displayed, enabling instructors to assess the full CPR cycle — compressions and ventilations — from a single unit without handling the manikin during the assessment. The standalone monitor format means this feedback capability can be added to the compatible W45117 manikin configuration without requiring a full manikin upgrade.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Gaumard Scientific products in India with over 30 years of experience supplying medical simulators and genuine OEM accessories to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

Frequently Asked Questions

Q: Which manikin is this CPR Monitor compatible with? A: The CPR Monitor (S180 / W45050) is compatible with Gaumard CPR manikin W45117. Confirm your manikin model before ordering.

Q: What parameters does this monitor display? A: Cadence (compression rate) and depth of cardiac compression, and airway ventilation — the complete CPR quality parameter set.

Q: How do I order and get support from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 to confirm compatibility before ordering.