



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



AED Trainer Plus 2

Item No. 1018143 [W46277]
Weight 3.15kg
Dimensions 35.6 x 30.5 x 30.5 cm
MPN: 8008-0050-01

[Read More](#)

SKU:

Categories:AED Trainers

Product Description

The AED Trainer 2 provides a variety of simulations to help responders become familiar with an AED device. This trainer provides simulated shock delivery and has no highvoltage capabilities, ensuring safety during training. It is pre-configured with 10 training scenarios that simulate realistic sudden cardiac arrest episodes. The trainer can also be used with any manikin. Software and manual are only available in English.

The **ZOLL AED Trainer Plus 2 (Item No. 1018143 / W46277)** is a professional-grade AED training device from ZOLL Medical, USA — a globally trusted name in resuscitation technology. Designed as a faithful replica of the ZOLL AED Plus, this trainer gives healthcare educators and first responder instructors the ability to deliver realistic, scenario-based AED and defibrillation training in a completely

safe, zero-voltage environment.

Key Features:

- **10 pre-programmed cardiac arrest scenarios** — realistic sudden cardiac arrest simulations covering a wide range of rhythms and clinical presentations
- **Simulated shock delivery** — replicates the full AED workflow including shock advisory and delivery, with no high-voltage output, ensuring complete safety during training
- **Voice prompt control** — instructor can control voice prompts to guide or challenge trainees through simulated rescue sequences
- **Compatible with any CPR manikin** — integrates seamlessly into existing BLS training setups without requiring brand-specific manikins
- **Authentic ZOLL AED Plus form factor** — trainees gain familiarity with the exact device they will use in real cardiac emergencies, reducing hesitation at the point of care
- **Durable construction** — suitable for sustained use in high-throughput training programmes
- **English-only software and manual** — designed for English-medium training environments

Ideal For:

- NMC/INC-compliant BLS and AED skills stations in medical colleges and nursing institutes
- ACLS, BLS, and first responder training programmes
- Hospital in-service training and code blue preparedness drills
- Emergency medicine, critical care nursing, and paramedic programmes
- Institutions standardising on ZOLL AED Plus devices for their clinical environments

Q: What device does the AED Trainer Plus 2 replicate? A: It is a training replica of the ZOLL AED Plus — one of the most widely deployed AEDs globally. Institutions that use ZOLL AED Plus devices in their clinical settings benefit especially from this trainer, as trainees gain direct familiarity with the exact device they will encounter in a real emergency.

Q: Is this device safe to use during training? A: Yes. The AED Trainer Plus 2 delivers only simulated shocks — there is no high-voltage output whatsoever. It is completely safe for use on trainees and manikins in a training environment.

Q: How many scenarios are available? A: It comes pre-configured with 10 training scenarios simulating a realistic range of sudden cardiac arrest episodes, providing instructors with varied presentations to work through with different trainee groups.

Q: Can it be used with any CPR manikin? A: Yes. The trainer is compatible with any CPR training manikin, so it integrates into your existing simulation lab setup without requiring additional investment in brand-specific manikins.

Q: Can the instructor control the scenario during training? A: Yes. The device gives the instructor the ability to control voice prompts during a simulated rescue, allowing them to guide, challenge, or adjust the pace of the training scenario in real time.

Q: Is this suitable for NMC-compliant skills labs in India? A: Yes. The ZOLL AED Trainer Plus 2 fulfils the AED operation and defibrillation training requirements for BLS/ACLS skills stations in NMC-compliant medical college skills labs. Its scenario variety also supports structured competency assessment across trainee cohorts.

Q: Is the device available in languages other than English? A: No. The software and user manual for this trainer are available in English only. Institutions requiring regional or Hindi-medium training prompts should consider alternative AED trainers in the SEM Trainers range.

Q: What is the warranty on this product? A: Warranty terms are as per ZOLL Medical standard policy. Contact SEM Trainers for details applicable to Indian institutional purchases.

Specification	Detail
Item No.	1018143 (W46277)
MPN	8008-0050-01
Weight	3.15 kg
Dimensions	35.6 × 30.5 × 30.5 cm
Make	ZOLL Medical, USA
Scenarios	10 pre-programmed
Language	English only
High Voltage	None – training use only