




SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Life/form® “Airway Larry” Adult Airway Management Trainer Head

Item No. LF03667U
Weight 0.99kg
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:BLS Adult

Product Description

Airway Larry® Airway Management Trainers simulate a non-anesthetized patient for practicing intubation, ventilation, suction and CPR techniques. They have realistic anatomy and landmarks including teeth, tongue, oral and nasal pharynx, larynx, epiglottis, arytenoids, false cords, true vocal cords, trachea, lungs, esophagus, and stomach. Choose one of those great value for money trainers

and present your students with the challenges they may face in the real world.

| | Head/torso | Ideal for | Digital intubation | Laryngospasm & Tongue edema | Stomach rise | Sucti-oning | Jaw thrust | Manual carotid pulse | E.T., E.O.A., P.T.L, L.M.A., Combitube®, and King System insertion | Additional features |
|--|------------------------------|--|--------------------|-----------------------------|--------------|-------------|------------|----------------------|--|--|
| Adult Airway Management Trainer Head <i>Item No:</i> 1012740 | male head | introductory classes | Yes | - | - | Yes | - | Yes | Yes | can be attached to any of the adult CRiSis® Manikins |
| Advanced Trainer Head <i>Item No:</i> 1020961 | male head mounted on a board | introductory and advanced classes | Yes | Yes | Yes | Yes | - | - | Yes | can be attached to any of the adult CRiSis® Manikins |
| Adult Airway Management Trainer <i>Item No:</i> 1005633 | male head mounted on a board | testing purposes and Combitube® practice | Yes | - | - | Yes | Yes | - | Yes | can be attached to any of the adult CRiSis® Manikins |
| Economy Adult Airway Management Trainer <i>Item No:</i> 1005780 | male head mounted on a board | rescue & advanced airway management training | - | - | - | Yes | - | - | only LMAs and Combitube® | durable construction and bifurcated lung |
| Adult Airway Management Trainer <i>Item No:</i> 1017871 | torso with a female head | difficult airway and CPR training | Yes | - | - | Yes | Yes | Yes | Yes | durable construction and bifurcated lung |
| Advanced Trainer Head <i>Item No:</i> 1019249 | male head mounted on a board | advanced airway management training | Yes | Yes | - | Yes | - | Yes | Yes | durable construction and bifurcated lung |

| | | | | | | | | | | |
|---|-------------------------------|--|-----|-----|-----|-----|---|-----|-----|---|
| S.A.L.A.D. Simulator <i>Item No:</i> 1021583 | male head | suction assisted laryngoscopy and airway decontamination | Yes | - | - | Yes | - | Yes | Yes | rigid suction catheter distracts lower mandible and tongue to permit laryngoscope insertion |
| Airway Management Trainer <i>Item No:</i> 1019928 | male torso | airway management and CPR training | Yes | - | - | Yes | - | Yes | Yes | durable construction and bifurcated lung |
| Adult Deluxe Airway Management Trainer <i>Item No:</i> 1018851 | male head mounted on a board | cricothyrotomy and advanced airway management | - | Yes | Yes | - | - | - | - | lifelike skin that will not be damaged by aggressive intubation |
| Adult Airway Management Trainer <i>Item No:</i> 1005774 | male torso mounted on a board | rescue airway and anatomy training | - | Yes | Yes | - | - | - | - | review of the pulmonary system and oral/nasal pharyngeal spaces |

The Life/form® "Airway Larry" Adult Airway Management Trainer Head (Item No. LF03667U | Also known as 1012740) is a compact, lightweight male head trainer weighing just 0.99 kg, manufactured by Nasco Healthcare, USA. Designed to simulate a non-anesthetized adult patient, it is an ideal entry-level airway management trainer for introductory classes in intubation, ventilation, suctioning, and CPR techniques.

The trainer features highly realistic anatomy and landmarks including teeth, tongue, oral and nasal pharynx, larynx, epiglottis, arytenoids, false cords, true vocal cords, trachea, lungs, esophagus, and stomach — providing trainees with a true-to-life experience of the challenges they may face in real clinical scenarios. It supports digital intubation, suctioning, manual carotid pulse palpation, and insertion of ET tubes, EOA, PTL, LMA, Combitube®, and King System devices. The trainer head can also be conveniently attached to any of the adult CRiSis® Manikins for integrated training. A downloadable user manual is available on the product page.

Q1. What is the Life/form® Airway Larry Adult Airway Management Trainer Head used for? It is used for practicing adult airway management techniques including intubation, suctioning, ventilation, manual carotid pulse palpation, and insertion of ET tubes, EOA, PTL, LMA, Combitube®, and King System devices on a realistic non-anesthetized patient simulation head.

Q2. What anatomical landmarks does this trainer feature? The trainer includes realistic teeth, tongue, oral and nasal pharynx, larynx, epiglottis, arytenoids, false cords, true vocal cords, trachea, lungs, esophagus, and stomach — providing an accurate and challenging training environment.

Q3. What airway devices can be practiced on this trainer? Trainees can practice insertion of endotracheal tubes (ETT), EOA, PTL, LMA, Combitube®, and King System airway devices, as well as digital intubation and suctioning.

Q4. Who is this trainer best suited for? This head trainer is ideal for introductory airway management classes in medical colleges, nursing institutes, EMS training programs, and ACLS/BLS certification centers looking for an affordable, lightweight, and realistic airway training head.

Q5. Can this head be attached to a full manikin? Yes, the Airway Larry Adult Airway Management Trainer Head can be attached to any of the adult Life/form® CRiSis® Manikins, enabling integrated airway and CPR training in a single setup.

Q6. How does this compare to other Airway Larry models? This entry-level head trainer supports digital intubation, suctioning, manual carotid pulse, and multi-device airway insertion, making it suitable for introductory classes. The Advanced Trainer Head (1020961) adds laryngospasm, tongue edema, and stomach rise for more advanced training, while the torso models add CPR and jaw thrust capabilities.

Q7. Is a user manual available? Yes, a downloadable user manual is available directly on the product page at semtrainers.com for easy reference and setup guidance.

Q8. What is the weight of this trainer? The trainer weighs just 0.99 kg, making it highly portable and easy to use in any classroom or simulation lab setting.

Q9. Where can I purchase this in India? You can enquire and order through SEM Trainers at **+91-88495 63724** or **+91-98791 03905**, or visit semtrainers.com.