



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Child Airway Management Trainers Replacement Stomach Only

Item No. LF03628U
Weight 0.011kg
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:Airway Management Child

Product Description

For use with Life/form® Child Airway Management Trainers

Item No. LF03628U (also W44138), 0.011 kg, Nasco Healthcare, USA. The companion consumable to LF03627U (done immediately prior). Stomach only — the lightest consumable in the child airway range at 0.011 kg. Keeping it tight and parallel to LF03627U format, with the direct cross-reference to the lung set.

DESCRIPTION (copy-paste into WordPress product description):

Why Choose the Life/form® Child Airway Management Trainers Replacement Stomach Only for Your Simulation Lab in India?

The Life/form® Child Airway Management Trainers Replacement Stomach Only (Item No. LF03628U / W44138) by Nasco Healthcare, USA is a 0.011 kg genuine OEM replacement stomach bladder for Life/form® Child Airway Management Trainers — the stomach-only consumable, separate from the replacement lung set (LF03627U). The inflatable stomach bladder in Life/form® Child Airway Management Trainers is the consequence-feedback mechanism that makes oesophageal intubation training clinically effective: when the tube is placed in the oesophagus rather than the trachea, the stomach inflates visibly while the lungs do not — providing the trainee with the immediate, unambiguous feedback that mirrors what happens in a real patient. Without an intact stomach bladder, this feedback is absent and the trainer cannot deliver the oesophageal intubation consequence training it is designed for.

The stomach bladder degrades independently of the lungs from repeated inflation cycles — developing micro-perforations that prevent consistent inflation, or losing the material compliance that produces visible distension. It must therefore be stocked and replaced on its own schedule. At 0.011 kg, this is the lightest consumable in the child airway range and can be stocked in multiples without logistics overhead. Only genuine OEM LF03628U material is guaranteed to fit the compatible child trainer airway modules with correct connection geometry and correct inflation characteristics. For simulation labs running regular paediatric intubation training sessions, maintaining a stock of LF03628U alongside LF03627U (replacement lungs) ensures the trainer can be restored to full training condition rapidly whenever either consumable degrades.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and genuine OEM consumables to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

Frequently Asked Questions

Q: Which Life/form® trainers is this replacement stomach compatible with? A: LF03628U is designed for Life/form® Child Airway Management Trainers by Nasco Healthcare. Contact us to confirm compatibility with your specific trainer configuration before ordering.

Q: Is the replacement lung set included in this kit? A: No. This is a stomach-only consumable. The replacement lung set for Life/form® Child Airway Management Trainers is a separate consumable (LF03627U) and must be ordered independently.

Q: How do I order and get support from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724. We supply genuine Nasco Healthcare consumables with pan-India delivery.