



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Infant Intubation Head

Item No. 101-115U
Weight 1.59kg
Dimensions 30 x 20 x 13 cm
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories: Airway Management Child

Product Description

This trainer of an infant allows practice and training of nasal and oral intubation. Features include: • Inflatable lungs and stomach • Vocal cords are highlighted in white for easy laryngoscope viewing • Anatomical landmarks such as uvula, vocal cords, glottis, epiglottis, larynx, arytenoid cartilage, trachea, esophagus.

Why Choose the Life/form® Infant Intubation Head for Your Simulation Lab in India?

The Life/form® Infant Intubation Head (Item No. 101-115U / 1017245) by Nasco Healthcare, USA is a 1.59 kg board-mounted infant intubation training head — 30 × 20 × 13 cm — and the only configuration in the Life/form® infant airway range that accepts 2.5 mm endotracheal tubes, the correct neonatal/small infant ET tube size. From the comparison table, this is the defining distinction from the Infant Airway Trainer with Stand (LF03623U) which accepts 4.0 mm tubes: the 101-115U is the correct

trainer when neonatal and small infant ET tube sizing practice is the training objective. Supplied with carrying bag and lubricant. A downloadable user manual is available on the product page.

The vocal cords are highlighted in white for easy laryngoscope viewing — the same design feature present on the Child Intubation Head (101-125U) — making correct blade placement and vocal cord identification visible to both trainee and instructor simultaneously during classroom laryngoscopy sessions. Anatomical landmarks include uvula, vocal cords, glottis, epiglottis, larynx, arytenoid cartilage, trachea, and oesophagus — the complete reference set for nasal and oral intubation technique training at infant proportions. Inflatable lungs provide visible chest rise confirmation of correct tracheal placement; the inflatable stomach bladder confirms oesophageal intubation with gastric distension. Suctioning technique training is supported. For simulation labs, NRP training programmes, neonatal nursing institutes, and NMC-mandated skills labs in India that specifically train on neonatal and small infant airway management with 2.5 mm tube sizing, this is the correct and only Life/form® platform at that ET tube specification.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators to 1,500+ institutions across 28 states. We provide pan-India delivery, on-site installation, faculty orientation, and dedicated after-sales support. All products carry genuine manufacturer warranty and meet ISO 9001 quality standards.

Frequently Asked Questions

Q: What ET tube size does this trainer accept? A: The Infant Intubation Head (101-115U) accepts 2.5 mm endotracheal tubes — the neonatal/small infant tube size. This is the key differentiator from the Infant Airway Trainer with Stand (LF03623U), which accepts 4.0 mm tubes.

Q: Does this trainer include a stand or carry case? A: No stand is included. The unit is supplied with a carrying bag and lubricant. A downloadable user manual is available on the product page.

Q: What is the warranty and how do I get support? A: All Nasco Healthcare products carry standard manufacturer warranty. For support, reach us at sem@semtrainers.com or call +91-88495 63724.

Supplied with carrying bag and lubricant.

	Head/torso	Age	Digital intubation	Tongue edema	Laryngospasm	Stomach rise	Suctioning	Manual carotid pulse	Defibrillation	Additional features
--	------------	-----	--------------------	--------------	--------------	--------------	------------	----------------------	----------------	---------------------

NG Tube Insertion Trainer <i>Item No:</i> 1017250	head mounted on a board	1-year old	-	-	-	Yes	-	-	-	trainer for instruction of nasogastric tube insertion and care
Infant Airway Simulator <i>Item No:</i> 1018863	full-body	1-year old	-	-	-	-	Yes	-	-	chest cavity containing realistic organs, nasal passage permits placement of NP tube
Infant Airway Trainer <i>Item No:</i> 1017954	head mounted on a board	1-year old	-	-	-	Yes	Yes	-	-	accepts 4.0 mm endotracheal tubes
Infant Intubation Head <i>Item No:</i> 1017245	head mounted on a board	1-year old	-	-	-	-	Yes	-	-	accepts 2.5 mm endotracheal tubes
Advanced Infant Intubation Head <i>Item No:</i> 1017236	head mounted on a board	1-year old	-	Yes	-	Yes	Yes	-	-	Vocal cords are highlighted for easy viewing
Airway Management Trainer <i>Item No:</i> 1014201	head mounted on a board	3-year old	-	-	-	-	Yes	-	-	accepts 4.0 mm endotracheal tubes
Advanced Child Airway Management Trainer <i>Item No:</i> 1019787	head	5-year old	-	Yes	Yes	-	Yes	Yes	-	can be attached to any of the Life/form® Child Airway Management Trainers or Child CRiSis Manikins

Child Airway Management Trainer Head <i>Item No: 1019788</i>	head	5-year old	-	-	-	Yes	Yes	Yes	-	can be attached to any of the Life/form® Child Airway Management Trainers or Child CRiSis Manikins
Advanced Child Airway Management Trainer <i>Item No: 1005698</i>	head mounted on a board	5-year old	Yes	Yes	Yes	Yes	Yes	Yes	-	also allows for E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and King System insertion
Basic Child CPR/Airway Management Torso <i>Item No: 1018865</i>	torso	5-year old	-	-	-	-	-	-	-	This torso is the foundation of the Child CRiSis™ system as arms, legs, and a defibrillation chest can be added later.
Advanced Child CPR/Airway Management Torso <i>Item No: 1018864</i>	torso	5-year old	Yes	Yes	Yes	Yes	Yes	Yes	Yes	also allows for E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and King System insertion