

**SEM****SEM TRAINERS & SYSTEMS**

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Catheterization Simulator, male

Item No. LF00855U
Weight 4.073kg
Dimensions 43 x 25 x 30.5 cm
Make: Nasco Healthcare, USA

[Read More](#)**SKU:****Categories:**Catheterization

Product Description

Catheterization is one of the basic skills all nursing students must learn, and to learn the techniques, they must practice. Any patient volunteers? Not likely! But the Life/form® simulator allows demonstration and practice in urinary catheterization without embarrassment or discomfort to the patient or the students. In addition, the catheterization simulator is useful for anatomical identification and demonstration of perineal care. With the re-designed Life/form® Male Catheterization Simulator, catheterization can be practiced just as on a real patient. A lubricated catheter can be inserted in the urethral orifice, passed through the urethra, and into the bladder. When the bladder is successfully entered, artificial urine (water) will flow from the catheter. The student will feel the normal restrictions caused by the mucosal folds, bulbous urethra, and the internal urethral sphincter, just prior to entrance into the bladder. The experience teaches proper positioning and movement of the penis to allow the catheter to pass easily with a minimum of discomfort to the patient. Just as in a real patient, the cuff

must be completely deflated before the catheter is removed. Improper use of a Foley catheter can result in damage to the simulator, just as when using a Foley catheter with a real patient. Includes a hard carry case, lubricant, silicone catheter, and instruction manual. Fluid supply stand not included.

Why Choose the Life/form® Catheterization Simulator, Male, for Your Nursing Skills Lab in India?

The Life/form® Male Catheterization Simulator (Item No. LF00855U, 4.073 kg, 43 × 25 × 30.5 cm, Make: Nasco Healthcare, USA) lets nursing students practice one of the most fundamental clinical skills — urinary catheterization — without the discomfort or embarrassment of practicing on a real patient. It's equally useful for anatomical identification and demonstrating perineal care.

With this re-designed simulator, a lubricated catheter can be inserted through the urethral orifice, passed through the urethra, and into the bladder, exactly as with a real patient. Students feel the normal anatomical restrictions — mucosal folds, the bulbous urethra, and the internal urethral sphincter — just before entering the bladder, teaching correct positioning and movement of the penis for smooth, minimal-discomfort catheter passage. When the bladder is successfully entered, artificial urine flows from the catheter, giving immediate confirmation of correct technique.

The simulator also enforces correct Foley catheter procedure: the cuff must be fully deflated before catheter removal, just as with a real patient, and improper technique can damage the simulator the same way it would cause injury in clinical practice — reinforcing correct habits through realistic consequence.

Included: Hard carry case, lubricant, silicone catheter, instruction manual.

Please note: a fluid supply stand is **not included**. To see the fluid-flow confirmation feedback described above, you'll need to connect your own fluid supply setup — contact SEM Trainers for recommendations if you need this sourced separately.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and clinical skills training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Q: Is the fluid supply stand included? A: No — it's sold separately. Without it, the simulator won't demonstrate the fluid-flow feedback on successful catheterization.

Q: What's included in the box? A: A hard carry case, lubricant, silicone catheter, and instruction manual.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.