



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724

Mobile 1: +91 98791 03905

Email: sem@semtrainers.com

Website: www.semtrainers.com



Catheterization Simulator PRO, female

Item No. 1020229 [P93P-F]

Weight 3.5kg


Dimensions 50 x 42 x 30 cm

Brand 3B Scientific

[Read More](#)

SKU:

Categories:Catheterization



Product Description

With the female P93 PRO catheterization trainer, female bladder catheterization procedures can be realistically demonstrated, practiced and assessed. The genital inserts are placed into the anatomically correct lower abdominal model and held in place with magnets. The material is soft and flexible: the labia are movable and can be separated. This way the student can practice all the necessary gestures

for catheterization (e.g. disinfection). The material gives a realistic feeling when inserting and removing the catheter into the bladder with realistic resistance and insertion depths. **Direct learning:** When the catheter has been correctly inserted, fluid runs out, just like with a real patient. The catheter can also be checked through the transparent bladder. The abdominal wall can be removed so that the process can be viewed better. **Includes suprapubic catheter:** If it is no longer possible to carry out catheterization through the urethra, a suprapubic catheter can be inserted. On the PRO version, a suprapubic catheter has already been put in place, and the student can practice cleaning it and taking care of it with the trainer. **All the benefits of the female P93 PRO catheterization simulator at a glance:**

- Female catheterization with realistic resistance
- Soft and movable labia (for practicing disinfection, for example)
- Anatomically realistic pelvic structure
- Liquid outflow if catheterization is successfully carried out
- The transparent bladder can be checked, and the abdominal wall can be removed
- Includes suprapubic catheter for practicing, cleaning and caring for it
- Stands securely on the table thanks to the non-slip feet
- Magnetic connectors for quick set-up and dismantling
- Easy to clean and maintain, and can be completely taken apart
- A sponge on the inside prevents a build-up of moisture

Delivery contents:

- Realistic lower abdominal model
- Catheterization set with transparent bladder
- Bottle and tube connectors
- Carry case

Note: For catheterization, we recommend using catheters with only one or two drainage eyes opposite each other. Both silicone and PVC catheters can be used. We recommend the following catheters:

- Urethral narrowing “open” = 16 CH catheter or 14 CH balloon catheter
- Urethral narrowing “partial narrowing” = 14 CH catheter or 12 CH balloon catheter

Catheters with a curved tip should not be used! Tip: According to the guidelines for correct catheterization, even during the simulation, lubricating gel (not included) should be placed into the urethra to ensure that the catheter can be inserted smoothly. Alternatively, some lubricating gel can be placed on the tip of the catheter.

Why Choose the Catheterization Simulator PRO, Female, for Your Skills Lab in India?

The Catheterization Simulator PRO, female (Item No. 1020229 [P93P-F], 3.5 kg, 50 × 42 × 30 cm, Brand: 3B Scientific) is a full anatomically realistic lower abdominal model for advanced female bladder catheterization training. Genital inserts mount magnetically into the abdominal model, with soft, movable, separable labia that let students practice the complete procedural sequence, including disinfection technique.

What distinguishes PRO from the BASIC trainer: the abdominal wall can be removed for a clearer view of catheter placement during teaching, and — uniquely — **a suprapubic catheter comes pre-installed**, so students can practice the ongoing care and cleaning of a suprapubic line, not just initial catheter insertion. This extends training beyond the urethral catheterization procedure into the downstream nursing care that follows when urethral catheterization isn't possible.

As with the BASIC model, successful urethral catheterization produces genuine fluid outflow, with placement visible through the transparent bladder. A sponge inside the model prevents moisture buildup, and the full unit can be taken apart for cleaning and maintenance.

Features:

- Female catheterization with realistic resistance
- Soft, movable, separable labia
- Anatomically realistic pelvic structure with removable abdominal wall
- Fluid outflow confirms successful catheterization
- **Pre-installed suprapubic catheter** for practicing care and cleaning
- Non-slip feet for secure table placement
- Magnetic connectors for quick setup/teardown
- Fully disassembles for cleaning; internal sponge prevents moisture buildup

Delivery contents: Realistic lower abdominal model, catheterization set with transparent bladder, bottle and tube connectors, carry case.

Recommended catheters (not included):

- Urethral narrowing "open": 16 CH catheter or 14 CH balloon catheter
- Urethral narrowing "partial narrowing": 14 CH catheter or 12 CH balloon catheter
- Use catheters with only one or two opposing drainage eyes. Both silicone and PVC catheters are suitable. **Do not use catheters with a curved tip.**
- Lubricating gel is recommended (not included) for smooth insertion.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of 3B Scientific products in India with over 30 years of experience supplying medical simulators and clinical skills

training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Q: How does PRO differ from BASIC? A: PRO offers a full anatomical abdominal model with a removable wall for better visualization, and comes with a suprapubic catheter pre-installed for care and cleaning practice — skills BASIC doesn't cover.

Q: Does this come with a carry case? A: Yes — the PRO version includes a carry case, unlike the BASIC model.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.