



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com



Life/form® Infant Male and Female Catheterization Trainer

Item No. LF01035U
Weight 0.75kg
Make: Nasco Healthcare, USA

[Read More](#)

SKU:

Categories:Catheterization

Product Description

The Life/form® Infant Male and Female Catheterization Trainer is an inexpensive 2-in-1 trainer that provides realistic practice of this difficult skills procedure. • Practice insertion, position, catheter balloon inflation and deflation • Feel resistance and pressure as with a real patient • Urine flows when catheter is in proper position • Practice both male and female infant catheterization training • Anatomically realistic with both male and female genitals • Genitals made of soft, lifelike material • Replaceable male and female genitals • Lower half of infant only • Design for use with a standard 8 French latex or silicone intermittent catheter Includes: infant manikin, one male and one female genital insert, 2 oz. lubricant, infant Foley catheter, and instruction manual.

Why Choose the Life/form® Infant Male and Female Catheterization Trainer for Your Pediatric

Nursing Skills Lab in India?

The Life/form® Infant Male and Female Catheterization Trainer (Item No. LF01035U, 0.75 kg, Make: Nasco Healthcare, USA) is an inexpensive, compact 2-in-1 trainer purpose-built for practicing this genuinely difficult pediatric procedure. Students practice catheter insertion, positioning, and balloon inflation/deflation, feeling realistic resistance and pressure throughout — with urine flow confirming correct catheter placement.

Swappable male and female genital inserts, made of soft, lifelike and replaceable material, let one compact lower-half infant manikin cover both training scenarios without needing two separate units. The trainer is designed for use with a standard 8 French latex or silicone intermittent catheter — a critical pediatric-specific sizing detail distinct from adult catheterization equipment.

Its small footprint and low cost make it a practical addition for pediatric nursing programs building out infant-specific clinical skills alongside adult catheterization trainers.

Includes: Infant manikin (lower half), one male and one female genital insert, 2oz lubricant, infant Foley catheter, instruction manual.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and pediatric clinical skills training equipment to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support.

Q: What catheter size does this use? A: A standard 8 French latex or silicone intermittent catheter.

Q: Does this cover both male and female infant catheterization? A: Yes — swappable soft genital inserts let one manikin cover both scenarios.

Q: How do I order from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724 for pricing and availability.